


# Using the Tefonius software



# USING THE Tafoneius SOFTWARE

Before the **Tafoneius** software can be used to control the ventilator, the ventilator needs to be initialised. This initialisation process is covered in the section “Getting Started with **Tafoneius**”

Once initialised, there are some common steps to take when preparing for a patient:

Pre-filling the system  
Entering patient information  
Choosing a Preset

<b>Starting a Session .....</b>	<b>6</b>
<b>The Finalise Recording procedure.....</b>	<b>7</b>
<b>Entering initial patient information .....</b>	<b>8</b>
<b>The Ventilator controls.....</b>	<b>12</b>
<i>Tidal Volume:</i> .....	13
<i>Respiratory Rate:</i> .....	13
<i>Inspiratory Time:</i> .....	13
<i>Inspiratory Pause:</i> .....	13
<i>Maximum Working Pressure Limit (MWPL)</i> .....	14
<i>CPAP/PEEP</i> .....	14
<b>Validating entries .....</b>	<b>14</b>
<b>Dependant Values: .....</b>	<b>15</b>
<i>The System button</i> .....	16
System Tab .....	16
System Errors .....	16
Windows Event Viewer .....	17
Leak/Compliance .....	17
Piston Initialisation .....	17
Leak & Compliance .....	17
<b>Running the Leak &amp; Compliance procedure .....</b>	<b>18</b>
<b>The Patient Data button.....</b>	<b>21</b>
<i>Data Logging</i> .....	21
<i>Finalise Recording</i> .....	21
<i>Data Storage Locations</i> .....	21
<b>Viewing Recorded Data .....</b>	<b>22</b>
<i>Viewing Recorded data during a recording</i> .....	22
Viewing the instantaneous value using Data View .....	22

# USING THE TAFONIUS SOFTWARE

<i>Viewing Recorded data from the Tafoneius Shell</i> .....	23
<b>Exporting Data</b> .....	<b>23</b>
<b>The Configure button</b> .....	<b>24</b>
<i>Presets</i> .....	24
<i>Initial Settings</i> .....	25
<i>User settings</i> .....	25
<i>Changing Presets</i> .....	26
Save As New .....	26
Save Active .....	26
Rename .....	26
Switch To .....	26
Delete .....	26
Reset Factory Default .....	26
<i>Gas Monitoring</i> .....	27
Gas monitoring tab with the IRMA Gas Analyser: .....	27
Gas monitoring tab with the Artema AION Gas Analyser: .....	28
<i>IRMA Setup</i> .....	29
<i>Pre-Use check</i> .....	31
<i>IRMA Monitor: General Description</i> .....	31
<i>Calibration, self-test and reset</i> .....	31
Oxygen Recalibration .....	31
CO <sub>2</sub> , Nitrous and Agent Calibration (Zero reference calibration) .....	32
<i>Using the ARTEMA AION Gas Analyser</i> .....	33
<i>Display of Measured Values</i> .....	33
<b>Other Control buttons: Alarm Silence &amp; Freeze Traces</b> .....	<b>35</b>
<i>Alarm Silence</i> .....	35
<i>Freeze Traces</i> .....	35
The countdown time for Freeze Traces is determined by the Countdown value set in the Initial settings attribute panel. ....	35
<b>Monitoring Software</b> .....	<b>36</b>
<i>Monitoring options</i> .....	36
<i>Organisation of traces</i> .....	37
Adding and Removing Traces .....	38
<i>Controlling the appearance of traces</i> .....	38
The ECG/IBP Attribute tab .....	40
The CO <sub>2</sub> Attribute tab .....	43
The Airway Pressure (AP) Attribute tab .....	44
The Oxygen Attribute tab .....	45


## USING THE Tafoneius SOFTWARE

The SPO2 Attribute tab.....	46
The Agent Attribute tab.....	47
The N2O Attribute tab.....	48
The Temperature Attribute tab.....	49
The Traces/Trends Attribute tab.....	50
<b>Alarms, Alerts and Critical Error Messages.....</b>	<b>51</b>
<i>Creating an Anaesthetic Record Chart.....</i>	<i>52</i>
Altering the Start Time of the Anaesthetic Record chart.....	52
<i>Viewing the Anaesthetic Record chart.....</i>	<i>52</i>
Data Recording using the Slow Data File.....	53
Changing the interval for Periodic Data in the Slow Data File.....	53
Effect of the Logging Interval on the Slow Data entries.....	53
<b>The Tafoneius Shell.....</b>	<b>56</b>
<i>Tafoneius.....</i>	<i>56</i>
<i>View Data.....</i>	<i>56</i>
<i>Export Data.....</i>	<i>58</i>
<i>Touch Screen Alignment.....</i>	<i>59</i>
<i>Recover Anaesthetic Records.....</i>	<i>59</i>
<i>Company Name.....</i>	<i>59</i>
<i>Shutdown.....</i>	<i>59</i>
<b>Appendix A : SLOW DATA INFORMATION.....</b>	<b>60</b>

# USING THE Tafonius SOFTWARE

## Pre-filling the system

After **Tafonius** has been initialised the piston is left at the bottom of the cylinder ready for filling. Make sure the Y-piece is capped, turn the vaporizer to the required pre-fill concentration and turn on the manual oxygen flow meter to fill the cylinder with gas. The cylinder icon on the left of the screen indicates the volume of gas in litres in the cylinder.



Fill the cylinder to a level equal to the maximum expected patient's tidal volume PLUS 2.5 Litres. e.g. 500kg horse. Average tidal volume predicted is 5.0L so a larger breath of 7.5L can reasonably be expected. Then add 2.5L, which is the buffer volume used by the ventilator. Total pre-fill volume = 10L. When the cylinder graphic shows it is at the required volume, turn off the oxygen flow.

This sets the volume of the "Virtual Bag" for spontaneous breathing. The concept of the Virtual Bag is central to how **Tafonius** operates. This action of the ventilator maintains a minimum system volume thereby reducing anaesthetic and oxygen consumption as well as decreasing the system's time constant.

### The Piston as a "Virtual bag"

When the ventilator is in standby mode and a patient is connected, the system behaves like a 'perfect' bag. When the patient breathes out there is no expiratory resistance and when the patient breathes in there is no inspiratory resistance. This ideal performance is accomplished with the Airway Servo System employed by Tafonius. The airway pressure is sensed at the centre of the Y-piece. As soon as a pressure deviation in excess of 0.5cm H<sub>2</sub>O from the ambient pressure is detected then the piston is moved so as to negate this pressure change. If the patient breathes out then the piston withdraws upwards, effectively filling "the bag". If the patient inspires then the piston moves down, effectively emptying. Because the piston is driven by a motor it is the motor that overcomes the resistances of any piping, valve or soda lime. To the patient it feels like breathing into room air albeit through and ET tube.

### Dictating the size of the "Virtual Bag"

With varying sizes of patients it is beneficial to keep the size of any reservoir bag in the system to a minimum. This then allows rapid responses to changes in anaesthetic or oxygen concentrations. The setting of the Tidal Volume and the concept of fixed buffer volumes control the size of the "Virtual Bag" in **Tafonius**.

tafonius employs two fixed buffer volumes; an upper and a lower. The lower buffer volume is 1.5L and the upper buffer volume is 1.0L. These volumes are constant and are added to the set Tidal Volume to determine the virtual "Bag" size thereby creating upper and lower volume limits.

The upper volume limit = Lower Buffer + Tidal Volume + Upper Buffer.

The lower volume limit = Lower Buffer

Whenever the piston volume exceeds the upper limit, gas is removed from the system. When the piston volume falls below the lower limit then gas is added to the system.

# USING THE Tafonius SOFTWARE

## Starting a Session

Tafonius uses the concept of anaesthetic or procedure **Sessions** for the collection and recording of data. As soon as the Tafonius software starts then information begins to be stored in the Slow Data file and the Session is started. This storage of data continues until a session is Finalised. A session may be Finalised by one of the following:

- Clicking or touching the Finalise Recording button in the Data Logging tab
- Exiting the program through the System and then Exit buttons
- Exiting the program by pressing and holding the Mains Control button
- An abrupt program closure such as power loss, closure of Windows without closing Tafonius (non-controlled shutdown)

In all but the last procedure, the data will be saved **and** associated with all the patient information.

In the event of a non-controlled shutdown there will be NO DATA LOSS, although the patient information will not be associated with the data. In this instance the Files associated with the session will be saved with the word "Patients" as a prefix and the Slow Data and Anaesthetic Record will use only the date as identifiers. For example, a session interrupted by a non-controlled shutdown on January 7<sup>th</sup> 2009 will have a record named "Patients2009Jan07" with Slow Data and Anaesthetic Records named AR2009Jan7-8\_50\_56 and SD2009Jan7-8\_50\_56 respectively. The suffix 8\_50\_56 indicates that the session began at 08:50:56 on the morning of January 7<sup>th</sup>.


For more information on entering Patient and case information, see the section entitled "Entering Initial Patient Information" later in this manual

For more information on the production of an Anaesthetic Record see the section entitled "Creating An Anaesthetic Record Chart" later in this manual

# USING THE TAFONIUS SOFTWARE

## The Finalise Recording procedure

Using any of the first three Finalising methods listed above will cause the following dialog to appear



**Accept Findings**

Do you accept these findings:

Hospital Number:	12345678
Patient Name:	Alfonzo
Owner Name:	Smithers
Patient ID:	192837465
Species:	Equine
Breed/Type:	Arab
Age:	4
Sex:	M
Body Weight:	500Kg
Notes:	Alfonzo notes

At Presentation:

Anaesthetist:	John Irving (House Surgeon)		
Surgeon:			
Procedure:			
Respiration Rate:	0 Breaths/min	Pulse:	0 Beats/min
Temperature:	0 Deg C	CRT:	1sec
PCV:		MM:	Pink
TTP:		Pant:	No
BUN:		Fasted:	No
Anaesthetic Risk Status:	1		

To accept these findings press Accept, to clear the finding press Discard, or to return to the application and change the findings press Change.

This dialog presents the current information associated with the recording session. If these are all correct then click or touch the Acept button. To edit any feature shown, click or touch the Change button. This will take you to back to the Patient Data edit region, where any of the details can be changed.


If all of the details are correct, then clicking or touching the Acept button will cause the following to happen automatically:

- The Session is ended.
- The patient details are used to create the appropriately named Slow Data and Anaesthetic Record files.
- A new Session is started using, either the next set of patient information in the Case list, or the default Patient/Anonymous details.
- Recording for the new Session is started

# USING THE Tafonius SOFTWARE

## Entering initial patient information

Touch the Patient Data button



A new window appears in the information area

A screenshot of the 'Patient Details' window. It has tabs at the top: 'Cases', 'Patient Details' (selected), 'Notes', 'At Presentation', 'Data Logging', and 'Data View'. The form contains fields for: 'Hospital Number:', 'Patient Name:' (with 'Anonymous' entered), 'Owner Name:', 'Patient ID (PID):' (with '111111' entered), 'Species:' (with 'Unknown' entered), 'Breed/Type:', 'Age:' (with '0' in 'y' and '0' in 'm' boxes), 'Sex:', and 'Weight:' (with '0' in 'Kg' box). A 'Close' button is at the bottom right.

Enter the Hospital Number (if used) and then use the tab key to move to the next data entry field. Fields that take numerical data entry will bring up the small on-screen numerical key-pad which can be used instead of a keyboard. By default the patient's name is "Anonymous". Edit this to enter your own patient's name. When the details have been entered touch the "At Presentation" tab.

A screenshot of the 'At Presentation' window. It has tabs at the top: 'Cases', 'Patient Details', 'Notes', 'At Presentation' (selected), 'Data Logging', and 'Data View'. The form contains fields for: 'Anaesthetist:', 'Surgeon:', 'Procedure:', 'Fasted' (radio buttons for Yes/No, with 'No' selected), 'Pant' (radio buttons for Yes/No, with 'No' selected), 'Anesthetic Risk Status' (radio buttons for 1, 2, 3, 4, 5, and E, with '1' selected), 'Temp' (radio buttons for Fahrenheit/Celcius, with 'Celcius' selected), 'RR:', 'PCV:', 'TP:', 'BUN:', 'Pulse:', 'CRT:' (with '1sec' selected), and 'MM:' (with 'Pink' selected). A 'Close' button is at the bottom right.

Fill in the procedure details along with the rest of the presenting clinical information.



## USING THE TAFONIUS SOFTWARE

Entries can be left blank. The Anaesthetist, Surgeon and Procedure fields are used when an Anaesthetic Record chart is produced. The Anaesthetist, Surgeon and Procedure entries are all taken from a Drop-Down list. To view the list click or touch the down arrow in the blue box next to the field you wish to enter details for. If the item you require appears on the list click or touch it to select it. This provides a quick way to enter commonly used information. If the item does not appear on the list then it may be added by touching or clicking on the ellipse button



Clicking or touching the Ellipse button opens an editor.  
Here the Anaesthetist Drop-Down list has been opened for editing:

Anaesthetists	
Dr Andrews	
Lorraine Gosch	
Dr Simpson	
John Irving (House Surgeon)	

Items in the list may be ordered using the Up & Down buttons. The Add, Edit and Delete buttons allow standard editing of the entries in the list. When editing is complete, click or touch the Close button to return to the “At Presentation” tab.

### Notes Tab:

Any notes that you want to enter relating to the case can be entered as free text in the Notes tab. These notes will appear in the Anaesthetic Record chart in the area under the Header of general information.

## USING THE Tafonius SOFTWARE

### Cases Tab:

Details relating to cases can be pre-entered as a list of patients. Associated case data can be entered and then recalled when required.

Patient Name	Patient ID (PID)
Anonymous	0000000
Alfonzo	192837465
→ Diamond Sunset	198725463

Select

New

Delete

By default, when the software is started details appear for a patient called Anonymous with a patient ID of 000000.

The green arrow next to the patient name shows that these patient details are in use. Patient details cannot be deleted if they are in use: you must first select another patient. Touching or clicking the New button opens the Patient details tab ready for new patient information. To change patients, click or touch the patient name or ID and then click or touch the Select button. The selected patient information is then used for the current recording session. When the monitoring/recording session is finalised, the Patient details in use at that time will be the ones used to create the recording files. During the Finalising procedure a dialog box will appear to check that the Patient details to be saved are correct.

The patient name and ID is always shown at the bottom of the Tafonius screen. Click or touch either of these two screen areas to open the Patient attribute dialog.

## USING THE Tafonius SOFTWARE

### Data Logging Tab:

In this tab the options for Slow Data Recording and Anaesthetic Record chart heart rate can be set. Because Tafonius has 3 means of monitoring heart rate (ECG, IBP & Pulse-Ox) it is necessary to select one to be used for the Anaesthetic Record chart. Unless you have specific reasons for changing the Slow Data Logging time, leave this set at 5 seconds.

The screenshot shows the 'Data Logging' tab selected in a software interface. The tab bar at the top includes 'Cases', 'Patient Details', 'Notes', 'At Presentation', 'Data Logging' (highlighted), and 'Data View'. The main content area contains the following elements:

- Slow Data Logging:** A text label followed by a text input field containing the number '5' and the unit 'sec'.
- Time:** A label in blue text.
- Start:** A text field showing '14:27:55'.
- Elapsed:** A text field showing '00:04:14'.
- Record Spontaneous Breaths:** A checkbox that is currently unchecked.
- Record Ventilated Breaths:** A checkbox that is currently checked.
- For Anaesthesia Record Log Heart Rate From:** A label in blue text above a group box containing three radio button options:
  - ECG:** The selected option, indicated by a green dot.
  - IBP:** An unselected option, indicated by a blue circle.
  - Pulse-ox:** An unselected option, indicated by a blue circle.
- Anaesthetic Record Start:** A button located at the bottom center.
- Finalise Recording:** A button located at the bottom right.
- Close:** A button located at the bottom right, outside the main content area.


When all the patient information has been entered, click or touch the "Close" button.

# USING THE TAFONIUS SOFTWARE

## The Ventilator controls

The ventilator controls consist of Tidal Volume, Respiratory Rate, Inspiratory Time, Maximum Working Pressure Limit and a CPAP (or PEEP) setting. These are the only controls that can be set for ventilation.

The main controls (Tidal Volume, Respiratory Rate & Inspiratory Time) have two sets of numbers associated with them. These numbers change depending on whether the system is in Standby (spontaneous breathing) mode or Ventilate mode. In Standby mode, the large blue numbers are the measured values derived from the patient's respiratory efforts. The Airway Servo System is working as a respiratory monitor. The smaller black numbers are the ventilator settings that will be used in Ventilate mode. If no breathing is detected after 60 seconds the blue values are all set to zero.



When you touch a button it changes colour and becomes ready to accept a new value. A numeric keypad appears in the lower right hand corner of the screen where you can enter a new value. Or the value may be changed by means of the UP/DOWN arrows. When the new value has been entered do one of the following:

- Touch the body of the button again,
- Touch the "Enter" button on the numeric keypad or
- Touch the body of the next button whose value you want to change.

## USING THE TAFONIUS SOFTWARE

### Tidal Volume:

In Standby mode: The set Tidal Volume is indicated by the smaller black number and is the volume that will be delivered per breath when the switch is made to Ventilate mode. It also sets the volume of the Virtual Bag used by the patient. The blue value indicates the breath-by-breath tidal volume of the spontaneously breathing patient

In Ventilate Mode: Sets the actual volume delivered to the patient. **Tafonius** uses automatic compliance compensation so the volume received by the patient is the volume entered. The volume displaced by the cylinder will be more than the set value to compensate for the volume lost to the compliance of the system. The patient actually receives the TV you set regardless of the length of the breathing systems tubes. This is very helpful when 20 or 30 ft tubes are connected to reach a patient in an MRI. Tidal Volume units are Litres

### Respiratory Rate:

In Standby mode: The blue figures show the measured respiratory rate, recalculated every 10 seconds on a rolling 1 minute basis. This gives a true reflection of respiratory rate. The smaller black figure shows the respiratory rate that will be used in Ventilate mode.

In Ventilate Mode: The large black number shows the respiratory rate in use.

Respiratory Rate units are Breaths per Minute

### Inspiratory Time:

The Inspiratory Time button controls the total length of the Inspiratory phase. There is also a smaller sub-button located within the Inspiratory Time button that allows a pause to be included in the Inspiratory time. This inspiratory pause is limited to 40% of the I time. It can be displayed in either seconds or percent of I time, by repeated touches of the sub-button.

In Standby Mode: The blue figures show the measured inspiratory time, recalculated after every breath i.e. on a breath to breath basis. The smaller black figure shows the total Inspiratory Time that will be used in Ventilate mode.

In Ventilate Mode: The large black figure indicates the total time for inspiration. This may not be the same as the Inspiratory flow time, depending on the setting for Inspiratory Pause.

### Inspiratory Pause:

In Standby Mode: Has no action

In Ventilate Mode: Controls the portion of the Inspiratory phase that is held as a pause. The maximum pause is limited to 40% of the total Inspiratory time. Increasing the Pause time will increase the time the inspiratory volume is held at the end of the Inspiration stroke. The total Inspiratory Time is unaffected so Inspiratory Flow rates will increase as more Pause time is added.

Touch the Pause button to activate it and open the numeric keypad. Touch it again to toggle between setting the pause as a % of Inspiratory time, or as a finite time. Enter the actual time or percentage required. Press "Enter" on the numeric keypad when done.

**In general, as with the ventilator settings above, any button shaded green on the display is a value you can set. These include the following.**

# USING THE Tefonius SOFTWARE

## Maximum Working Pressure Limit (MWPL)

The Maximum Working Pressure Limit control sets a level, above which ventilation is prevented. This is particularly useful when dealing with neonates or during open-chest surgery when lung expansion is not restricted by a closed thorax. The Maximum Working Pressure Limit is an absolute value and is unaffected by any settings of PEEP or CPAP. If the MWPL value is reached during an Inspiratory phase then the ventilator immediately sounds an alarm and changes to the expiratory phase, where once again the cylinder behaves as a Virtual Bag. The MWPL setting can be useful for patients that are “fighting” the ventilator to prevent excessive pressures developing.

## CPAP/PEEP

In Standby Mode: CPAP or Continuous Positive Airway Pressure can be applied during spontaneous breathing. With CPAP patients inspire and expire normally with no added effort, although the end expiratory pressure is elevated to the CPAP setting. This aids in maintaining open alveoli. The airway pressure is held constant during all phases of breathing.

In Ventilator Mode: PEEP or Positive End Expiratory Pressure applies during IPPV. Ventilating pressures rise during the Inspiratory phase as normal, but are allowed only to fall as far as the PEEP setting during expiration

## Validating entries

All entries are validated once complete. This means that any invalid entries that would exceed minimum or maximum settings that **Tefonius** is capable of delivering will not be accepted. When this occurs a warning sound is issued and a description of the input error is shown at the bottom of the screen as a black on yellow text message. In most circumstances the limiting factor preventing the input is also shown. Here, for example a value of 1.0 second has been entered as a Pause Time during an Inspiration Time of 2.0 seconds. This represents 50% of the Inspiratory Time, thus exceeding the maximum of 40% and so is disallowed. The maximum allowable entry of 0.8 seconds is used.

Increasing the *I<sub>t</sub>* to 2.5 seconds would allow a 0.5 second pause. The Pause Time box is outlined in red to indicate the source of the error.

Tidal Volume 5 0 (Litres)	Respiratory Rate 5 0 (Breaths per Minute)	Inspiratory Time 2 0 (seconds)	Pause 0.8 s
Minute Volume 25 0 (litres)	Inspiratory Flow 250 (litres/minute)	I:E Ratio 1:5 1:0	Expiratory Time 10 0 (seconds)

The Inspiratory Pause Time setting 1s was out of range. Valid range is 0-0.8s. Change limited by IT.

## USING THE Tafonius SOFTWARE

### **Dependant Values:**

The four buttons below the ventilator controls show the derived values for Minute Volume, Inspiratory Flow, I:E Ratio and Expiratory Time. They are dependent upon the settings of TV, RR and Itime


In Standby Mode: The large blue figures show the values calculated from the measurements made during normal inspiration and expiration. Respiratory rate is calculated on a rolling minute average that is updated every 10 seconds. The measured minute volume is an integral of the ventilated volume over the last minute. Again this is updated every 10 seconds. If no respiration is detected during a period of 60 seconds then these blue values are all reset to zero. The smaller black figures show the calculated values based on the values set for the three independent ventilator settings (TV, RR & IT).

In Ventilate Mode: The large black figures show the calculated values based on the values set for the 3 independent controls (TV, RR & IT).

# USING THE Tafonius SOFTWARE


## Control buttons: Controlling the Software

There are three buttons that control the look, feel and behaviour of the **Tafonius** software. These three buttons – System, Patient Data & Configure are located at the bottom right hand side of the screen.



### The System button

The System button opens a two-tab screen that contains a System and a Configuration tab. The Configuration tab can be entered directly from the “Configure” button on the lower right hand side of the screen.



### System Tab

The System tab contains buttons for System Errors, Leak/Compliance test settings, Language and User Manual options. Currently the Language is set to English only. Other language options will become available in due course. The User Manual button also presents options which at the time of writing are not yet available.

### System Errors

Any error or fault condition that occurs in the **Tafonius** software is reported in this System Errors list. An error is indicated by the flashing of the border of the System button in red and an audible alarm. The purpose of the list is to indicate the source of an unexpected alarm or fault condition and to assist with trouble-shooting. Some errors are transient and will self-clear, in which case no action is



## USING THE Tafonius SOFTWARE

necessary. Some errors will not self-clear or will remain for example after a device has been removed. In these instances, select the error and then click or touch on the “Resolve” button. If the error is no longer present then it is resolved and disappears from the list.

### Windows Event Viewer

This is reserved for use as a trouble-shooting tool. A field service engineer or someone from technical support may ask you to use this button to view any error sequences that have arisen during the use of Tafonius. Do not use this button during normal use as there may be an effect on the performance and behaviour of the Tafonius software.


### Leak/Compliance

A leak test and Compliance measurement is offered whenever the **Tafonius** software is launched. However sometimes it is convenient or necessary to run either utility without having to close the program and start again. Clicking on the Leak/Compliance button starts a sequence of events that proceeds in the following order:

(If the piston is not initialised then Piston Initialisation must be done before either leak testing or compliance measurement.)

### Piston Initialisation

This involves making sure that there is nothing attached to the Y-piece before the piston is sent to the bottom of the cylinder to identify and reset its zero position. If the piston has not been initialised since turn-on then the cylinder representation on the screen will have a hatched background. Once the piston has been initialised this changes to a solid green background. Leak testing and compliance values can then be measured.



Uninitialised piston background is hatched. Image also has a red flashing border

Initialised piston background is green and shows piston volume in litres




### Leak & Compliance

This button begins the process of first leak testing and then measuring the compliance of the system. For leak testing the piston is taken to the mid-point of the cylinder and then the piston produces a pressure of 20cm H<sub>2</sub>O. During this pressure phase the piston movement and hence leak is measured. The measured leak should be below 1500mls/min and in normal circumstances is less than 500 mls/minute. Values above or around 1500 mls/minute represent a serious leak and should be investigated. However, due to the way **Tafonius** maintains pressures and controls piston movement, even leaks of 1500mls/minute are unlikely to affect ventilator performance. They do however represent a loss of gas and agent.


# USING THE TAFONIUS SOFTWARE

## Running the Leak & Compliance procedure


Press the Leak/Compliance button on the System tab to begin the procedure




A dialog appears announcing that it is establishing a connection with the Auxiliary controller. This can take 5-10 seconds. After that a second dialog appears asking that any patient be removed from the circuit. When the breathing circuit is ready, press the Leak & Compliance button



At this point a stopper must be applied to the Y-piece to seal the system. Then press the Leak & Compliance button.




The piston is moved to the mid-point and then stops briefly before establishing a system pressure of 20cm H2O.




## USING THE Tefonius SOFTWARE

If for any reason there is an excessive leak the following message appears:




In this instance seal the leak and press the Retry button.

If the leak test is successful the leak value is displayed and the system then proceeds to perform the two compliance measurements.



A compliance measurement is made at the top cylinder position and at the 5.0L cylinder position and uses the leak value obtained in the first test. At the end of the test the leak and both compliance values are shown. Note that the values obtained will vary markedly depending on the nature of the breathing circuit connected. The two compliance values are used to obtain the equation for the straight-line graph of compliance for your machine and the breathing system in use at that time. This equation is then used to compute the compliance compensation required at any piston position and pressure.

## USING THE Tefonius SOFTWARE



Press "Finish" to leave the Leak & Compliance test procedure.

# USING THE **Tafonius** SOFTWARE

## **The Patient Data button**

Entering patient data has been described under the section “Entering Initial Patient Information”. This covers the “Patient Details” and “At Presentation” tabs found once the “Patient Data” button is pressed. Also found here are the “Data Logging” and “Data View” tabs.

### **Data Logging**

All physiological data and many automatic events are recorded in a file called the Slow Data (SD) file. How often data is entered in this SD file is dictated by the setting in the Data Logging tab. The default value is 5 seconds. A range between 1 second and several minutes can be entered.

Some natural events such as inspiration and expiration can be entered into the SD file. This however can fill the file with information that may not necessarily be required. Spontaneous breaths or ventilated breaths can be selected or deselected for entry in the SD file by checking or un-checking the respective tick-boxes.

The time that a particular recording session began and its duration are also shown in this tab.

### **Finalise Recording**

As soon as the **Tafonius** software starts, data is recorded in the SD file. When a patient session is ended the recording should be finalised. To do this push on the “Finalise Recording” button. A dialog box may appear asking to save the current Preset (see the section on User Presets) and then a further dialog box will appear stating that the data will be committed and that no further changes can be made. Click or touch OK to commit the data and begin a new session. When a session is Finalised, those patient details are removed from the Patient List in the Cases tab.

### **Data Storage Locations**

Where the data is stored will depend whether any Patient Information has been entered.

If no Patient data has been entered then the data is stored in a folder that is named using the word “Patients”, the calendar month and the day. Several patients or sessions may be held in this folder if no patient information is entered during any of the day’s operations.

e.g. An operation on September 12<sup>th</sup> 2008 where no patient information was entered will be stored in a folder called “Patients2008Sep12”

If patient data has been entered then the folder name will be derived from the Patient’s name and the patients ID.

e.g. Wild Billy , ID 1786394 will have a session folder created called WildBilly\_1786394. If there is more than one session recorded for Wild Billy on that day then all sessions will be stored in this folder.

Recorded data can be easily viewed and or exported using the utilities supplied with the **Tafonius** software. See the sections “Viewing Recorded Data” and “Exporting Data”


# USING THE TAFONIUS SOFTWARE

## Viewing Recorded Data


Recorded data can be viewed by two methods: During a recording session in the ventilator control software or after a session in a spreadsheet viewer.

### Viewing Recorded data during a recording

To view data on-the-fly, click or touch on “Patient Data” and then the “Data View” tab.




A trend viewer pops up showing the trend of the selected parameter from the beginning of the recording session to the current time. This data is constantly updated at a rate set by the Slow Data logging period. Use the drop-down list to select a different parameter.



### Viewing the instantaneous value using Data Drop Down List

Click or touch a point in the graph area of the Data View tab. A line appears at the point of touch. At the bottom of the line on the x-axis of the graph is the time at which that point of data was recorded. Along the right hand side of the line are the instantaneous value(s) of the associated data.



# USING THE **Tafonius** SOFTWARE

## **Viewing Recorded data from the Tafonius Shell**

Click or touch on the “View Data” button in the **Tafonius** Shell screen. This opens a Record Picker box to choose the Data Record to view. Data is presented in date order with the most recent at the top. Pick a record to view and then click/touch the “Select” button. A file or list of files will be shown that begin with SD. These are the Slow Data files recorded during a session. Touch or click or touch on the one required and then press the “Select” button. A spreadsheet viewer opens to display the recorded data. Touching a cell within a column representing measured physiological data will display that data automatically as a trend graph.

To create graphs for other parameters that are not physiological such as the set Tidal Volume or Compliance Compensated Tidal Volume it is necessary to export the data to a spreadsheet viewer of your choice. This is easily done from the **Tafonius** Shell.

## **Exporting Data**

From the **Tafonius** Shell, click or touch on the “Export Data” button. An “Export Data” dialog box opens listing the recorded data on the **Tafonius** computer. Data is presented in date order with the most recent at the top. Choose the Record folder required and then click or touch on the “Next” button. Then choose the destination for your exported data. Typically this will be a USB flash drive inserted in the hub at the front of the screen. When the destination has been set, click or touch on the “Export” button to copy the data across to the selected device. Data is exported as a Comma Separated Value (CSV) file. The Anaesthetic Record chart in pdf format will be exported at the same time.

# USING THE Tafonius SOFTWARE

## The Configure button


The “Configure” button is used to access the following: Presets, Initial Settings, Gas Monitoring and SPO2 setup. The “Security” button in the Configure window is not active at this time.

### Presets

Presets are a means of individually saving **Tafonius** software layouts and initial settings that suit a particular purpose. For example you may have a preset called “Foals” or a preset called “Colics” or even “Senior Anaesthetist”.

Each Preset saves the way the screen is configured including how the traces are arranged as well as all the initial values used when the program starts, as well as all the alarm settings and whether they are enabled or not.

There are two ways in which Presets are used to load values at the start of a procedure. There are the Initial Settings, which describe the preferred initial values for such items as Tidal Volume, Respiratory rate etc and there are the user settings which remember the preferred settings for e.g. ECG speed and trace position. The values for Initial Settings can be changed by clicking on the Configure button at the bottom right hand side of the screen. A new smaller window area appears below the traces area.



Choose the preset you want to change by clicking on it in the Presets list. If it is already in use it will have the word Active after the preset name. If not, then click or touch on the required Preset and then choose the “Switch To” button. The selected Preset will be loaded. Then click or touch on the Initial Settings button on the left hand side.



# USING THE Tafonius SOFTWARE

## Initial Settings

Using the keyboard or the on-screen keypad enter the **initial** values you would like in each of the named fields.

The fields have the following meanings:

**TV:** Tidal Volume in litres.

**RR:** Respiratory rate in breaths per minute

**It:** Total Inspiratory time

**MWPL:** Maximum Working Pressure Limit

**IP:** Inspiratory Pause in seconds

Show IP as % ☐ : Shows the Inspiratory pause as a percentage of the **It**.

**CPAP:** Continuous Positive Airway Pressure in cm H<sub>2</sub>O

**FiO2:** Desired initial oxygen fraction. (not available in software version 1.0)

**APNOEA:** Length of time allowed between breaths before Apnoea alarm sounds

**Countdown:** Length of countdown timer before pop-up menus close automatically, in seconds

**Countdown** ☐ : Enable the countdown timer for closing pop up menus

**Mouse** ☐ : Enable the mouse cursor on the screen

When done, click or touch on the Preset tab again and then choose "Save Active". The next time this Preset is loaded these new initial values will be used.

## User settings

During the course of using the **Tafonius** program you may make changes to the screen appearance or to the appearance of a certain monitor trace. When changes such as these are made, that affect the personal way that the screen is viewed the Preset will be marked as altered and an asterisk will appear next to the name in the title bar.



At any time or at the end of a procedure this unique set up of Tafonius can be saved by again choosing to save the Preset. If this is not done explicitly during use then you will be prompted to save the Preset when finalising a patient session. You can choose to change the active preset by touching SAVE or you may save the current setup as a new setup by touching SAVE AS NEW and renaming it as desired.


Note that you will not be prompted to save a changed Preset when you shut down the Program using the main ON/OFF switch.

# USING THE Tafoneius SOFTWARE

## Changing Presets

When the **Tafoneius** software starts it uses, by default, the last Preset in use. To change a Preset, click or touch on "Configure". The Preset list is displayed. Along the bottom of the Preset Window are 5 buttons, used to control how Presets are created, named, selected and deleted plus two other buttons to Reset Factory Defaults or Close the Configuration dialog.

The Preset currently in use is indicated by the word Active in brackets after the Preset name. This name is also displayed at the top of the screen in the title bar.



### Save As New

This button creates a copy of the Active Preset and names it "Custom Preset". Use the Rename button to change the name to one of your choice. If you have made changes to the stored active preset the new stored preset will include those changes.

### Save Active

This button saves any changes made to the Active Preset

### Rename

This button allows a Preset to be renamed. Type the new name in the edit box that appears and click or touch "OK".

### Switch To

Use this button to change to a different Preset from the Active one. To use this button, click or touch on the Preset you wish to use and then click or touch on the "Switch To" button.

### Delete

This button removes the highlighted Preset. Note that it is not possible to delete an active preset.

### Reset Factory Default

This button will apply all the default factory settings to the Active Preset. Note that the Preset must be Active for the changes to be made to it. A confirmation dialog box appears stating which Preset will be changed. If this is the Preset that you wish to Reset, then click or touch "OK".


# USING THE TAFONIUS SOFTWARE

## Gas Monitoring

This tab displays all of the status and diagnostic information for the current Gas Monitor in use. At the current time a mainstream device called IRMA or a sidestream device called Artema AION is used for CO<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>O and Agent monitoring. Your machine will be fitted with one or other of these units.

### Gas monitoring tab with the IRMA Gas Analyser:

Clicking on the Gas Monitoring tab reveals a Window with 5 further tab options: Status, Sensor, Errors, Config & Data



#### Status tab

Used for fault diagnostic information only

#### Sensor tab

Used for fault diagnostic information only

#### Errors tab

Used for fault diagnostic information only

#### Config tab

Used for fault diagnostic information only

#### Data tab

Location of the Zero option for the IRMA monitor. See the section on "Use of the IRMA Monitor" for full information on when and how to use the "Zero" option found in this tab.

In most situations it is only necessary for users to use the Data tab. Information from the other tabs may be requested by a service engineer or a technical support person during fault diagnosis.

## USING THE TAFONIUS SOFTWARE

Gas monitoring tab with the Artema AION Gas Analyser:

The screenshot shows the Tafonius software interface. On the left is a vertical sidebar with buttons: 'Presets', 'Initial Settings', 'Gas Monitoring' (highlighted in yellow), 'SpO2', and 'Security'. The main window has a top bar with 'Configuration' and 'System' tabs. Below this is a sub-menu with 'Modes/Status' (highlighted), 'Status', 'Alarms', 'Cmds', 'Demo', and 'Calib'. The 'Modes/Status' panel contains a 'Reset' button, 'Op Mode' set to 'ISOAccuracy', 'Comm Mode' set to 'Normal', 'Water Trap' set to 'Adult', 'Atm (hPa)' set to '102.1', 'Pump' set to 'Regulating', and 'Flow' set to '202'. A 'Message' field at the bottom displays 'Command OK'. A 'Close' button is located in the bottom right corner of the main panel.

### Modes/Status

Used for fault diagnostic information only

### Status

Used for fault diagnostic information only

### Alarms

Used for fault diagnostic information only

### Cmds

Used for fault diagnostic information only

### Demo

Used for demonstration purposes only.

### Calib

Used for fault diagnostic information only

In most situations it is only necessary for users to use the Cmds tab for performing a zero calibration. Information from the other tabs may be requested by a service engineer or a technical support person during fault diagnosis.

## USING THE Tafonius SOFTWARE

### Using the PHASEIN IRMA Monitor

#### Important Note:

*The IRMA mainstream multi-gas probe is intended to be connected to a patient breathing circuit for monitoring of inspired/expired gases of patients in intensive care, anaesthesia and emergency care. It is not intended to be used as the only means of monitoring a patient. It shall always be used in combination with other vital signs monitoring devices and/or professional human judgements of patient condition.*

*The IRMA multi-gas analyser is intended to be used by trained and authorised veterinary professionals only. It is not intended to be used in outdoor transport applications such as in cars or in aircrafts.*



#### IRMA Setup

**Plug the IRMA device into the IRMA socket located on the bulkhead behind and below the utility drawer, on the Right Hand Side of the machine.**

**Snap the IRMA sensor head on top of the IRMA airway adaptor. It will click into place when properly seated**



**A green LED indicates that the IRMA sensor is ready to use**



## USING THE Tefonius SOFTWARE

Connect the 15mm male IRMA/Airway adaptor connector to the breathing circuit. When used with a foal 22mm circuit connect the IRMA device as shown in the next 2 images.



Connect the 15mm female IRMA/Airway adaptor connector to the patient's endotracheal tube



When used with a 50mm mainstream adaptor the IRMA adaptor is integrated into the Large airway adaptor as shown:



## USING THE TAFONIUS SOFTWARE

**Always position the IRMA sensor with the O2 cell (or top) pointing upwards**



### Pre-Use check

Before connecting the IRMA adaptor to the breathing circuit, verify the O2 calibration by checking that the O2 reading on the monitor is correct (21%). See the following section “**Oxygen Recalibration**” on how to perform room air calibration.

### IRMA Monitor: General Description

The IRMA breath monitor unit is a small cube-shaped unit housing a replaceable Airway Adaptor and Oxygen sensor. The main housing for the unit has an LED in one corner which indicates the status of the sensor. These states are described in detail in the section "IRMA Status & Error Messages". The IRMA unit connects to the **Tafonius** monitoring system via a 9-Way Female D Type connector. When the unit is connected and the unit powered ON, the IRMA device runs through a self-test sequence. This is indicated by the LED on the main housing changing alternately through GREEN, RED, BLUE a number of times. Once the initialisation sequence is complete the LED stays GREEN. If the IRMA unit has the capability of detecting anaesthetic agents then the LED will turn BLUE in the presence of these agents.

IRMA units may be connected/disconnected while powered without any problems.

### Calibration, self-test and reset

The IRMA unit is self-calibrating and does not require scheduled calibration on a monthly or even annual basis. Calibration is achieved in a number of ways but must be done while the unit is plugged in and powered ON.

**ALL CALIBRATION MUST BE DONE WITH THE UNIT IN FREE AIR AND NOT WHILST CONNECTED IN A PATIENT CIRCUIT**

### Oxygen Recalibration

For calibration of the Oxygen sensor, removal of the Airway adaptor forces an automatic recalibration to 21%. For this reason it is important that the removal and refitting of the airway adaptor is done in conditions where the airway concentration is at 21% such as in a well ventilated room.

1. Remove the Airway Adaptor by pulling the white plastic connector away from the main housing.
2. Wait until the LED flashes RED and the alarm sounds on the monitor.
3. Replace the Airway Adaptor.
4. Check that the LED turns green
5. Check that the O2 reading on the Monitor is 21%

## USING THE TAFONIUS SOFTWARE

### CO<sub>2</sub>, Nitrous and Agent Calibration (Zero reference calibration)

These can all be reset from the Windows software via the Configure/Gas Monitoring/Data tab route. Only perform this function when asked to do so by the Windows program or when indicated by the following conditions:


- 1) A repeated "Gas Concentration out of Range" error message that appears when measuring room air and is not cleared by disconnecting and re-connecting the IRMA device.
- 2) Incorrectly reported high Inspiratory CO<sub>2</sub> concentrations

**ALL CALIBRATION MUST BE DONE WITH THE UNIT IN FREE AIR AND NOT WHILST CONNECTED IN A PATIENT CIRCUIT**

**THE IRMA UNIT MUST BE ALLOWED TO WARM UP AND BE RUNNING FOR AT LEAST 15 MINUTES BEFORE A ZERO REFERENCE CALIBRATION IS PERFORMED**

To RESET the IRMA device, click or touch the "Zero" button in the Data tab of the Gas Monitoring Window found in the Configuration Setup.

A warning dialog box will appear with the following message:



Make sure the conditions for re-calibration are correct before continuing, then click or touch Yes. If you are unsure about the conditions required for calibration, click or touch "No" and read the above paragraphs carefully before running the "Zero" option again.



## USING THE TAFONIUS SOFTWARE

### Using the ARTEMA AION Gas Analyser

The Artema AION Gas Analyser is a sidestream gas analyser, drawing the sampled gas from the Y-piece. Other than the physical differences and the inherent latency between breathing and waveform display, the reporting and display of all gas parameters is the same as with the IRMA Gas Analyser.



Sidestream sampling port (capped)

The sidestream sampling line connects to the sampling catheter that enters the Y-piece as shown. The sampling catheter extends down into the ET tube for optimum sampling.

### Display of Measured Values

THE Gas Analyser (IRMA or AION) will continuously display and update the numerical and trace information for the patient. The trace region shows the capnogram as either a line or filled trace. The main numerical data region associated with the trace shows the inspired and expired values of CO<sub>2</sub>.



The larger figure represents the expired value and the smaller figure shows the inspired value. The third value to the right is the respiratory rate in breaths per minute. If a period of apnoea occurs then the alarm will sound and the respiratory rate will be replaced by a red box, which displays the number of seconds elapsed since the last breath.



## USING THE Tafoneius SOFTWARE

Clicking or touching this red box will silence the Apnoea alarm. The box will remain red and continue to count the seconds since the last recorded breath. A fresh breath will clear the Apnoea state and return to displaying Respiratory Rate (RR).

During periods of Apnoea, the inspired values are constantly updated indicating the gas concentrations in the breathing circuit.

DRAFT

# USING THE TAFONIUS SOFTWARE

## Other Control buttons: Alarm Silence & Freeze Traces

### Alarm Silence

Press the button once for a 5-second silence of all alarms. After 5 seconds the alarms, if active, will again sound.



Press and hold the Alarm Silence button until it shows the legend "Locked".



The audio portion of all current alarms are then suppressed. Any new alarm condition will override the locked out condition and again sound the alarm.

### Freeze Traces

Use this button to freeze all trace information on the screen. All other monitoring activities continue as normal. The button will change colour to indicate that it is active.



Also a countdown timer will begin showing on the button. At the end of the countdown period the freeze option will be lifted. Alternatively simply clicking or touching the Freeze button once more will release the freeze option.

The countdown time for Freeze Traces is determined by the Countdown value set in the Initial settings attribute panel.

# USING THE TAFONIUS SOFTWARE

## Monitoring Software

### Monitoring options

Tafonius has the options of monitoring the following parameters via its in-built software:

Single Channel ECG

Single Channel Invasive (Direct) Blood Pressure

Breath Analysis:

Inspired/Expired CO<sub>2</sub>

Inspired/Expired O<sub>2</sub>

Inspired/Expired N<sub>2</sub>O

Inspired/Expired Anaesthetic Agent:

Halothane

Isoflurane

Sevoflurane

Desflurane

Enflurane

Single Channel Pulse-Oximeter


Dual Channel Temperature with delta

Airway Pressure

# USING THE Tafonius SOFTWARE

## Organisation of traces

The monitoring screen consists of 3 main regions – a trace region where all time-traces are displayed, a numerical region where all values are displayed and a status region where status information about **Tafonius** is displayed.




The trace region is centralised on the screen and in the picture above contains the ECG, CO2, Airway Pressure and Pulse-Oximetry trace. This area can be organised as required by placing any trace or combination of traces here. The trace region itself is divided into three horizontal areas each of which is one trace high. It is possible to place a trace in any of these 3 horizontal areas. Additionally it is possible to have multiple traces in any one of these horizontal areas.

## USING THE Tefonius SOFTWARE

### Adding and Removing Traces

To change the trace screen viewing options touch any trace. A new editing box appears at the bottom of the screen




The text in the yellow strip indicates the trace that has been touched (selected). Here the user has two options: Touching inside the Trace/Trend dialog box will remove the trace from the screen and place it in the reserve area. Alternatively touching the screen in another area of the trace region will place that trace alongside and to the left of any other traces in that horizontal area.

Any monitor that has been placed in the reserve area can be put back on the screen in any of the three horizontal areas at any time. To do this touch any "number box" e.g. the ECG heart rate box. The set of tabbed attribute panels appear at the bottom of the screen. Select the Traces/Trends tab and then select and drop the desired monitor trace back on the screen. "Pick and put" the desired trace where you want it.

### Controlling the appearance of traces

All traces have various attributes associated with them that determine for example their height, speed, shift and alarm options amongst others.

To access these options for any monitor touch the number panel associated with the monitor. For example to access the options for the ECG monitor, touch the left hand side of the ECG/IBP panel at the top right of the screen.



## USING THE TAFONIUS SOFTWARE


The monitoring attribute panels appear at the bottom of the screen opened at the ECG tab. Note that the ECG attribute tab shares its place with the IBP attribute tab as they are both located on this common ECG/IBP tab. When the left hand side of the ECG/IBP panel is touched, the ECG attributes are shown. If the right hand side of the ECG panel is touched, which has the IBP values, then the attribute tab will be opened with the IBP attributes showing. To change between ECG and IBP attributes touch the appropriate radio button at the left hand side of the ECG/IBP tab

All alarm settings and trace settings are particular to any given Preset and so the Preset should be saved after making adjustments to any of the controls in the monitoring attribute panels.

## USING THE TAFONIUS SOFTWARE

### The ECG/IBP Attribute tab

There are two active sections to this tab: Traces and Alarms. In the Traces section all the controls relate to how the trace can be viewed. Values can be set for either ECG or IBP options by selecting the appropriate option in the Select box.



### Trace options: ECG

**Sensitivity:** controls the height of the ECG trace seen on the screen


**Sweep:** Controls the screen sweep speed in rates of mm/second.

**Filter:** Allows the setting of a digital filter to remove unwanted signal noise. Use a notch filter for mains interference. Use a Low-pass filter for tremors due to movement e.g. shivering. Move the slider in the Filter box to control the degree of filtering.

**Shift:** Shifts the base-line of the ECG up and down the screen. Use the shift option to centralise the ECG for optimum viewing.

**Heart Rate:** Select the box next to "Beep on ECG" to produce an audible beep with every detected ECG complex. Since **Tafonius** only has one audio output device, only one monitor can be made to beep at any time. Selecting ECG for example, will mean that any previous selection is overridden.

To change between ECG and IBP attributes touch the appropriate radio button at the left hand side of the ECG/IBP tab



### Synchronising the ECG and IBP waveforms

ECG and IBP waveforms can be synchronised in time as long as they have the same trace width on the screen. Normally they would be placed one above the other so that a visual association between the



## USING THE Tafonius SOFTWARE

two traces can be made. If the two traces are not synchronised then select the “Synchronise” option in the ECG/IBP tab. If the two traces do not have the same width they cannot be synchronised and the option will be greyed.

### Trace options: IBP

The screenshot shows the 'ECG/IBP' tab selected. The 'Trace' button is highlighted. The 'Alarms' button is below it. The 'Select' section has radio buttons for 'ECG' and 'IBP', with 'IBP' selected. The 'Synchronise' checkbox is checked. The 'Sensitivity' section has a vertical slider with values 2, 5, 10, 20, 30, and 60 mmHg/div. The 'Sweep' section has radio buttons for 12.5, 25, 50, 100, and 200 mm/sec. The 'Shift' section has a horizontal slider. The 'Heart Rate' section has a 'Beep On IBP' checkbox. The 'Zero IBP' button is below the 'Synchronise' checkbox. The 'Close' button is in the bottom right corner.

### Sensitivity

Moving this slider will change the full scale height of the IBP trace in accordance with the legend on the right hand side of this scale. Note that there is an automatic levelling feature in the IBP trace that means that the bottom of the IBP trace is adjusted after every sweep to lie on the base line. This will always keep the IBP trace in view even if the blood pressure values are slowly changing.

### Sweep

Selecting one of the radio buttons will set the sweep speed for the IBP trace. Note that this option may be further controlled or changed depending on the state of the Synchronise check box and, or the presence of ECG and IBP traces on the Tafonius screen. If the ECG and IBP signals are not synchronised then the sweep speed of ECG and IBP traces can be different. If the traces are synchronised then the IBP sweep speed will be set to and follow the ECG sweep speed.

### Shift

The Shift slider alters the position of the baseline of the IBP trace.

### Heart Rate

Select the box next to “Beep on IBP” to produce an audible beep with every detected pulse waveform. Since **Tafonius** only has one audio output device, only one monitor can be made to beep at any time. Selecting Beep On IBP for example, will mean that any previous beep selection is overridden.

### Zero IBP:

When a direct arterial sampling line has been placed and the sensor positioned at the level of the patient’s heart (but with connection to ambient pressure not arterial pressure) the “Zero IBP” button is used to set the current output to zero.

## USING THE TAFONIUS SOFTWARE

### Alarm options: ECG & IBP

Click or touch the Alarms button on the ECG/IBP tab to show the alarm options for ECG and IBP monitors.

Because the monitor is intended for use with only one animal at a time there are potentially several sources of heart rate data, not least of which are ECG and IBP. When setting a heart rate alarm it is possible to choose to derive the heart rate from either ECG or IBP depending which is in use or which is more reliable in the current circumstances.

The screenshot shows the 'Alarms' settings screen for the ECG/IBP tab. The screen has a top navigation bar with tabs: ECG/IBP, CO2, Oxygen, Agent, N2O, SpO2, AP, Temp, and Traces/Trends. The 'ECG/IBP' tab is selected. On the left, there are two buttons: 'Trace' and 'Alarms', with 'Alarms' being highlighted. Below these buttons is a checkbox labeled 'Synchronise' which is checked. The main area is titled 'Alarm Limits' and contains four sections: 'Heart Rate', 'MAP', 'Systolic', and 'Diastolic'. Each section has 'High' and 'Low' value fields, an 'Enable' checkbox, and a 'Source' selection (ECG or IBP). The 'Heart Rate' section shows a High value of 55 bpm and a Low value of 0 bpm, with the ECG source selected. The 'MAP', 'Systolic', and 'Diastolic' sections all show High and Low values of 0 mmHg, with the Enable checkboxes unchecked. A 'Close' button is located in the bottom right corner.

### Alarm Limits

This group of alarms includes alarms for Heart Rate, Systolic, Diastolic and Mean Arterial Pressure. To set any alarm value, click or touch the associated edit box. The numerical keypad will appear on the right of the screen. Use this or the main keyboard to enter the numerical value and then click or touch Enter.

### Heart Rate

Upper and lower alarm values can be set for the heart rate and these can be independently enabled by touching or clicking the Enable box next to the appropriate Alarm edit box.

### Systolic

Set high and low limits for the systolic blood pressure value. High and Low alarms can be independently enabled.

### Diastolic

Set high and low limits for the diastolic blood pressure value. High and Low alarms can be independently enabled.

### Mean Arterial Pressure

Set high and low limits for the mean arterial blood pressure value. High and Low alarms can be independently enabled.

# USING THE TAFONIUS SOFTWARE

## The CO2 Attribute tab

### Properties

**Units:** Set the displayed units as Volume %, millimetres of mercury (mmHg) or as kilopascals (kPa). When units are changed, there is an automatic change in the units applied to alarms so that the corresponding value in the new units is used.

**Trace:** Select whether to display the Capnogram as a filled trace or as a simple line curve.

**Scale:** Fixes the vertical scale of the CO2 trace.

**Beep on Breath:** Select this tick box to produce an audible beep with every detected breath.

### Alarm Options


To set Alarms associated with CO2 measurement, click or touch the Alarms button

				Enable
Inspired	High	---	%	<input type="checkbox"/>
Expired	High	---	%	<input type="checkbox"/>
	Low	---	%	<input type="checkbox"/>
Breathing Rate	High	---	bpm	<input type="checkbox"/>
	Low	---	bpm	<input type="checkbox"/>

To set an alarm, place a tick in the enable box and then enter the appropriate alarm value.

## USING THE TAFONIUS SOFTWARE

### The Airway Pressure (AP) Attribute tab



#### Sensitivity


This slider sets the scales of the Airway Pressure trace. The value next to the pointer represents the maximum positive height of the trace window.

#### Sweep

This slider sets the sweep speed of the Airway Pressure trace, which is the speed at which the trace is drawn across the screen.

#### Airway Pressure Alarms

To access the Airway Pressure Alarms, click or touch the Alarms button



The CPAP alarm allows an alert condition when the measured CPAP pressure deviates outside of the required CPAP value. Alarm units are in percentage values of the set CPAP value.

## USING THE TAFONIUS SOFTWARE


### The Oxygen Attribute tab

ECG/IBP CO2 **Oxygen** Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Sweep**

Slow  Fast

**Scale**

☐ 0-50%

☒ 0-100%

Close

#### Sweep

Use this slider control to set the sweep speed of the Oxygen trace across the screen.

#### Scale

Use this control to set the full scale height of the Oxygen trace window. Note that there is no shift feature with this trace. To view Oxygen values in excess of 50%, the 0-100% scale range must be used.

#### Alarm Options

Click or touch the Alarms button to show the Alarm options for Oxygen.

ECG/IBP CO2 **Oxygen** Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Alarm Limits**

			Enable
Inspired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
Expired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
FiO2	+/-	<input type="text"/> %	<input type="checkbox"/>

Close

The Oxygen Alarm options allow the alarm limits for Inspired and Expired oxygen to be set. Alarms can be set for High and Low Inspired and Expired values. Click or touch the Enable box to enable the alarm and then enter the alarm value by clicking or touching the appropriate alarm edit box.

**FiO2** - This option is not yet implemented and will be available in the future.

# USING THE TAFONIUS SOFTWARE

## The SPO2 Attribute tab

ECG/IBP CO2 Oxygen Agent N2O **SpO2** AP Temp Traces/Trends

Trace

Alarms

Sweep

Slow Fast

Sound

☐ Beep On SpO2 Pulse

Close

### Sweep

Use the Sweep slider control to change the speed at which the SpO2 trace moves across the screen

### Sound

Click or touch the box next to “*Beep On SpO2 Pulse*” to enable audible beeps associated with each detected pulse waveform. This will over-ride any other beep selection.

### Alarm Options

Click or touch the Alarms button to show the Alarm options for SpO2.

ECG/IBP CO2 Oxygen Agent N2O **SpO2** AP Temp Traces/Trends

Trace

Alarms

Alarm Limits

				Enable
Pulse Rate	High	<input type="text"/>	bpm	<input type="checkbox"/>
	Low	<input type="text"/>	bpm	<input type="checkbox"/>
Saturation	Low	<input type="text"/>	%	<input type="checkbox"/>

Sound

☐ Beep On SpO2 Pulse

Close

### Alarm Limits

Click or touch the enable box to enable the alarm and then enter the appropriate value.

### Sound

Click or touch the tick box to select audible beeps every time a pulse wave on the SPO2 monitor is detected.

## USING THE TAFONIUS SOFTWARE

### The Agent Attribute tab

ECG/IBP CO2 Oxygen **Agent** N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**In Use**

- ☒ No Agent
- ☐ Halothane
- ☐ Enflurane
- ☐ Isoflurane
- ☐ Sevoflurane
- ☐ Desflurane

**Sweep**

Slow  Fast

**Scale**

- ☒ 0-5%
- ☐ 0-10%

Close

**In Use:** The IRMA OR detector used for agent monitoring does not have automatic agent ID so it is necessary to set the agent being used. Select the appropriate radio button next to the agent in use. This change will be reflected in the Agent button on the screen.

### Sweep

Use the Sweep slider control to change the speed at which the Agent trace moves across the screen

### Scale

Use this control to set the full scale height of the Agent trace window. Note that there is no shift feature with this trace. To view Agent values in excess of 5%, the 0-10% scale range must be used.

### Alarm Options

Click or touch the Alarms button to show the Alarm options for Agent.

ECG/IBP CO2 Oxygen **Agent** N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Isoflurane Alarm Limits**

			Enable
Inspired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
Expired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>

Close

**Alarm Limits:** Set the alarm limits required for a particular agent here. Note that different alarm limits can be set for each and every agent. Click or touch the Enable box to enable the alarm and then enter the appropriate value.

## USING THE TAFONIUS SOFTWARE

### The N2O Attribute tab

ECG/IBP CO2 Oxygen Agent **N2O** SpO2 AP Temp Traces/Trends

Trace

Alarms

Sweep

Slow Fast

Scale

☐ 0-50%

☒ 0-100%

Close

### Alarm Options

Click or touch the Alarms button to show the Alarm options for N2O.

ECG/IBP CO2 Oxygen Agent **N2O** SpO2 AP Temp Traces/Trends

Trace

Alarms

Alarm Limits

Inspired	High	<input type="text" value="---"/>	%	Enable <input type="checkbox"/>
Expired	High	<input type="text" value="---"/>	%	Enable <input type="checkbox"/>

Close

### Alarm Limits:

This tab allows the alarm limits to be set. Alarms can be set for High Inspired and Expired values. Click or touch the Enable box to enable the alarm and then enter the alarm value by clicking or touching the appropriate alarm edit box.



## USING THE TAFONIUS SOFTWARE

### The Temperature Attribute tab

ECG/IBP CO2 Oxygen Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

Time

- ☐ 10min
- ☐ 20min
- ☒ 1hr
- ☐ 2hrs
- ☐ 5hrs
- ☐ 10hrs

Temperature Scales

- ☐ Fahrenheit
- ☒ Celcius

Close

#### Time

Selecting an option here sets the amount of data displayed in the trace window. By selecting a long Time period, trends can easily be visualised. By choosing a short time period, details in trace changes are more easily observed.

**Temperature Scales:** Choose to use either Fahrenheit or Celsius as the temperature scale. Alarm values change accordingly.

#### Alarm Options

Click or touch the Alarms button to show the Alarm options for Temperature.

ECG/IBP CO2 Oxygen Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

T1 T2 Delta

Alarm Limits

		Enable
High	<input type="text"/>	<input type="checkbox"/>
Low	<input type="text"/>	<input type="checkbox"/>

Alarm Silence


Close

**Alarm Limits:** There are 3 alarm options: Alarms for T1 temperature probe, Alarms for T2 temperature probe and Alarms for the Delta temperature. (The Delta temperature is the difference in temperature between the two temperature readings of T1 and T2). Click or touch the enable box to enable the alarm and then enter the appropriate value.

## USING THE TAFONIUS SOFTWARE

**Alarm – Silence:** Use this button to silence a temperature alarm. An alarm will sound when a temperature probe is discovered to be absent. When starting the software if no temperature sensors are present an alarm may sound and the border of the temperature gauge will flash red. Touch the temperature gauge to open this Attribute tab and use the Silence button to silence the alarm.

### The Traces/Trends Attribute tab



For a detailed description of how traces are managed and moved around the screen, see the section entitled “Adding and Removing Traces” at the start of the section above.

## **USING THE Tafoneius SOFTWARE**

### **Alarms, Alerts and Critical Error Messages**

DRAFT

# USING THE Tafonius SOFTWARE

## Creating an Anaesthetic Record Chart

The Tafonius software creates an automatic record that can be viewed as an Anaesthetic Record Chart. Every 5 minutes monitoring data is recorded and entered on the chart in a similar manner to a manual record. Information about the patient, surgeon, anaesthetist and other details can all be recorded.


All non-clinical information is entered via the Patient Data button. This process is covered fully in the section entitled "Entering Initial Patient Information"

During a patient Session, there are effectively two recorded sets of data. One set of data is saved as the Slow Data file. This records periodically (default is every 5 seconds) all of the information relating to the running of Tafonius. For a more details explanation of the Slow Data files see the section entitled "Data Recording using the Slow Data File".

The other set of data is data saved as the Anaesthetic Record file. This file consists of data saved every 5 minutes as plotted or tabled data as well as header information such as name of Anaesthetist or Surgeon.

## Altering the Start Time of the Anaesthetic Record chart

It is not uncommon to prepare Tafonius some time in advance of connection to the patient. During this time there will be continual recording of non-meaningful data. To avoid this, the Start time for the Anaesthetic Record data can be set using the Anaesthetic Record Start button. Pressing this button will reset the start time for the Anaesthetic Record data ONLY. Slow data information is not affected. The Anaesthetic Record Start button is found in the Data Logging tab of the Patient Data dialogs. The following dialog will appear:



This procedure can be repeated as many times as required, noting that each time, data prior to selecting this option will no longer be available.

After this creation of the Anaesthetic Record chart is automatic and the chart can be viewed from the Tafonius Shell via the View Data button.

## Viewing the Anaesthetic Record chart

The Anaesthetic Record chart can only be viewed from the Tafonius Shell, whilst Tafonius is **NOT** running. Click or touch the View Data button in the Tafonius Shell. Use the Record Picker to select a patient and then click Select. Then highlight the AR file and click or touch Select. This will open a pdf version of the Anaesthetic Record in Adobe Viewer.

# USING THE TAFONIUS SOFTWARE

## Data Recording using the Slow Data File

The Slow Data file is an automatically generated file that records all information gathered during a Session. Data contained in this file includes all the monitoring information in terms of numerical and trend data. This periodic data (Default is every 5 seconds) also contains all events that occur during use. If an alarm condition occurs this is recorded in the Slow Data file. Similarly a change of anaesthetic agent type or the level of the MWPL setting is recorded in the Slow Data file. This raw file is non-editable so that the user can be sure that all events that occurred are recorded. For a full description of the data saved in the Slow Data file, see Appendix A, Slow Data Information.

The Slow Data file can be exported to a Flash disk using the Export data function accessed from the Tefonius shell. See the section entitled "Export Data" later in this manual.

In addition to automatic events, user events can be added to the Slow Data file. These can be added by clicking or touching the Events button at the bottom right hand side of the Tefonius screen. These events will appear in the Slow Data file as well as being used for the events shown in the Anaesthetic Record chart. For more information on entering manual Events see the section entitled "Adding Events During a Procedure" later in this manual.

## Changing the interval for Periodic Data in the Slow Data File

The interval for periodic data can be changed in the Data Logging tab of the Patient Data dialog. Click or touch the Patient Data button and then select the Data Logging tab. In the Slow Data Logging edit box enter the value in seconds to be used as the logging interval. Note that this value will have an effect on the final file size. Logging every second will create a substantially larger file than every 5 seconds. Note also the effect of the logging interval as explained below.

## Effect of the Logging Interval on the Slow Data entries


Whatever the interval period between logging of entries, automatic event entry will always be logged. For example, if the Ventilate button is pressed an automatic event will be generated and added to a temporary list of events. When the next logging episode occurs, all events stored since the previous logging event will be entered into the SD file. This is different from the physiological monitoring data which is simply grabbed at the instant of the periodic logging.

Care should be used when interpreting periodic logged data because events that appear at the same logging could be separated by up the periodic logging interval.

## USING THE TAFONIUS SOFTWARE

### Adding Events during a procedure

Events can be added to the Slow Data (SD) record during a procedure. Events can be pre-existing in a drop down list or may be added as a new text entry. To add an event, click or touch on the “Events” button at the bottom right hand side of the screen.



This will open a new editor panel at the bottom of the screen.

A screenshot of the 'Event List' editor panel. It contains a table with three columns: 'Date/Time', 'Event Description', and 'Status'. The first three rows are populated with data: '06/01/2009 15:50 - Fluid Administration:', '06/01/2009 15:49 - Intubation:', and '06/01/2009 15:49 - Anaesthetic Induction:'. There are 'Add' and 'Close' buttons on the right side.

Date/Time	Event Description	Status
06/01/2009 15:50	Fluid Administration:	
06/01/2009 15:49	Intubation:	
06/01/2009 15:49	Anaesthetic Induction:	

Click or touch on the “Add” button and then click or touch on the arrow of the Standard Events drop down box. In the example below, Anaesthetic Induction has been chosen from the Drop Down List.

A screenshot of the 'Event Details' editor panel. It shows a date and time selector set to '05 November, 2008' and '17:40:58'. Below this is a text input field containing 'Anaesthetic Induction:'. There are 'Accept', 'Clear', 'Standard Events:', 'Anaesthetic Induction' (in a dropdown menu), 'Add/Remove Standard Events', and 'Cancel' buttons.

Event Details: 05 November, 2008 17:40:58

Anaesthetic Induction: |

Standard Events: Anaesthetic Induction

Clear Add/Remove Standard Events Cancel

## USING THE TAFONIUS SOFTWARE

Standard Events can be added or removed from the Drop Down List by using the “Add/Remove Standard Events” button. Clicking or touching this button opens a new editor which lists all the current Standard Events. These can be Edited or Deleted using the Edit/Delete buttons. To add new items to the Standard Event list, click or touch the Add button and then type your new Event title.

The screenshot shows a window titled "Standard Events". On the left, there are two buttons: "Up" and "Down". In the center is a list of standard events: "Start Procedure", "Fluid Administration", "Premed", "Anaesthetic Induction", "End Anaesthesia", "Arrival At Recovery Site", "Extubation", and "Sternal". To the right of the list are four buttons: "Add", "Edit", "Delete", and "Close".

Once a new Event has been entered it will appear last on the list. To change the position of this new item in the list, use the UP and Down buttons to move the item. Once the new Event and the order of Events has been set, click or touch the Close button to leave the Standard Events editor.

To add your own entry click or touch in the Event Details box and start typing

The screenshot shows a window titled "Event Details:". At the top, there are two dropdown menus for date and time, showing "05 November, 2008" and "17:50:45". Below these is a large text input area with the placeholder text "Add a custom event here...". To the right of this area is an "Accept" button. At the bottom left is a "Clear" button. To its right is a "Standard Events:" label followed by a dropdown menu. Below this dropdown is a button labeled "Add/Remove Standard Events". To the right of the "Standard Events:" section is a "Cancel" button.

Click or touch “Accept” to add it to the SD file.

Use the “Clear” button to remove an entry from the Event Details edit box.

# USING THE Tafonius SOFTWARE

## The Tafonius Shell

The Tafonius Shell is the starting point for using the Tafonius Ventilator software and also the starting point for various related utilities.



The buttons:

### **Tafonius**


Launches the Tafonius ventilation and monitoring software

### **View Data**

This opens a browser for selecting a patient. All data records are shown in date order with the most recent data at the top of the list. Choose the Patient Record you would like to View and click or touch the Select button. Each patient will have a Slow Data file (SD file) and an Anaesthetic Record file (AR file).



## USING THE TAFONIUS SOFTWARE



To view the Slow Data information, highlight the appropriate file and then touch or click the Select button. The Slow Data is viewed in a spreadsheet format and can be used to show monitoring trends. By touching any of the measured parameter columns (e.g. ECG heart rate, Pulse-Ox saturation) a trend view of the data over the course of the procedure is shown at the bottom of the screen. Touching the column again hides the trend view. For greater versatility, export the data to a USB stick and view the data in a commercial spreadsheet of your choice. All Slow Data files are saved as CSV or comma separated value files.


To view the Anaesthetic Record chart, highlight the desired AR file and then touch or click Select. A new window opens to show the Anaesthetic Record pdf file.

For printing and archiving of either the Slow Data or the Anaesthetic Record use the Export Data button from the Tafonius Shell to save the data to a USB flash drive and then open the files in your host computer.

## USING THE TAFONIUS SOFTWARE

### Export Data

This button opens a browser to allow you to choose the data to export. All data in the browser is listed in date order with the most recent files at the top of the browser.




First select the patient that you wish to export data for and then click or touch the Next button. Then choose your destination drive (usually a USB flash drive named D:\) and touch or click the Export button. The Slow Data and the Anaesthetic Record data will be copied to your chosen destination.

## USING THE Tafonius SOFTWARE

### Options

This opens a small window with 3 facilities.



### Touch Screen Alignment

Touching or clicking this button starts a utility to align the touch screen to your computers screen. Simply touch the red circles in the way you would normally use the touch screen and follow the prompts to reset the alignment. If you make a mistake simply wait and the utility will time out and be cancelled. You can then try again.

### Recover Anaesthetic Records

Touching or clicking this button will create a new Anaesthetic Record chart if the original has been destroyed or corrupted for some reason.

### Company Name

Enter your practice or company name here. This will appear on the Tafonius shell as well as at the top of the Anaesthetic Record chart.

### Shutdown

This will shut down the PC and put the PC in an off condition. The PC can be started again by either turning on the Tafonius ventilator or by pressing the large ON/OFF button on the side of the PC.

# USING THE TAFONIUS SOFTWARE

## Appendix A : SLOW DATA INFORMATION

Monitoring data recorded in the Slow Data File includes the following:

### User control settings

CPAP/PEEP	Maximum PEEP or CPAP reading cm H <sub>2</sub> O
TV	Set Tidal Volume L
TVcc	Compliance compensated tidal volume L
RR	Respiratory Rate – breaths per minute
IT	Inspiratory Time in milliseconds
IP	Inspiratory Pause time in milliseconds

### Resultant parameters calculated from user settings

MV	Minute Volume L
I Flow	Inspiratory Flow L/min
I:E Ratio	Inspiratory:Expiratory ratio
Exp Time	Expiratory time in milliseconds

### Measured during Standby

TV	Spontaneous breath Tidal Volume in mls
Insp Time	Inspiratory time in ms
Exp Time	Expiratory time in ms

### Resultant parameters calculated from MEASURED values during Standby

RR	Respiratory rate in breaths per minute
MV	Minute Volume in litres
Iflow	Inspiratory flow rate in L/min
I:E Ratio	Inspiratory:Expiratory ratio

### Airway Pressure Measurements

CPAP/PEEP	CPAP(Standby) or PEEP(Ventilate) values in cm H <sub>2</sub> O
PIP	Peak Inspiratory Pressure in cm H <sub>2</sub> O

### Gas analysis Measurements

Insp CO <sub>2</sub>	Fractional %Volume measurement of inspired CO <sub>2</sub>
Exp CO <sub>2</sub>	Fractional %Volume measurement of expired CO <sub>2</sub>
Resp Rate	Breathing rate in breaths per minute
Insp Agent	Inspired Agent Fractional %Volume measurement
Exp Agent	Expired Agent Fractional %Volume measurement
Insp O <sub>2</sub>	Inspired Oxygen %Volume measurement
Exp O <sub>2</sub>	Expired Oxygen %Volume measurement
Atmospheric Pressure	Current Atmospheric pressure in mm Hg

### Nonin SpO<sub>2</sub> Measurements

SPO <sub>2</sub> HR	Heart rate from pulse-oximeter device
Saturation	Percentage oxygen saturation of arterial blood

## USING THE TAFONIUS SOFTWARE

### IBP Measurements

Mean	Mean arterial blood pressure (mean 4 average)
Systolic	Systolic blood pressure (mean 4 average)
Diastolic	Diastolic blood pressure (mean 4 average)
Heart Rate	Patient heart rate from blood pressure device

### Temperatures

T1	Temperature from Thermometer 1 in either °C or °F
T2	Temperature from Thermometer 2 in either °C or °F

### ECG HR

Heart rate from ECG monitoring device

### Battery Data

Batt1	Voltage level of battery 1 (lower 12v battery)
Current1	Charging rate of battery 1 in amps
Batt2	Voltage level of battery 2 (upper 12v battery)
Current2	Charging rate of battery 2 in amps

### Piston Position


Units	Shows the piston position in absolute units
-------	---

In addition to these headings in the file, there will be additional events recorded. The table below shows some typical events that have been recorded.

[08:09:33.21 - Preset Loaded: Custom Preset]  
[08:09:33.191 - Leak 600.0, Compliance Low 27.4, Compliance High 63.3]  
[08:09:33.201 - MWPL: 35]  
[08:09:33.311 - Cylinder Open]  
[08:09:33.702 - FMSO ON: 65535ms]  
[08:09:33.722 - Init Requested]  
[08:09:35.695 - Auxiliary Controller - No Watchdog]  
[08:09:35.705 - Auxiliary Controller - No Position data]  
[08:09:35.895 - Auxiliary Controller - Watchdog data established]  
[08:09:37.237 - FMSO ON: 65535ms]  
[08:09:37.247 - Init Requested]  
[08:09:37.778 - Auxiliary Controller - Position data established]  
[Thursday, March 27, 2008 08:09:38] Records Started - TAFONIUS Version 0.1.1092.0.

Every entry has a time associated with it and describes the event that occurred.

# Using the Tefonius software



# USING THE Tafoneius SOFTWARE

Before the **Tafoneius** software can be used to control the ventilator, the ventilator needs to be initialised. This initialisation process is covered in the section “Getting Started with **Tafoneius**”

Once initialised, there are some common steps to take when preparing for a patient:

Pre-filling the system  
Entering patient information  
Choosing a Preset

<b>Starting a Session .....</b>	<b>6</b>
<b>The Finalise Recording procedure.....</b>	<b>7</b>
<b>Entering initial patient information .....</b>	<b>8</b>
<b>The Ventilator controls.....</b>	<b>12</b>
<i>Tidal Volume:</i> .....	13
<i>Respiratory Rate:</i> .....	13
<i>Inspiratory Time:</i> .....	13
<i>Inspiratory Pause:</i> .....	13
<i>Maximum Working Pressure Limit (MWPL)</i> .....	14
<i>CPAP/PEEP</i> .....	14
<b>Validating entries .....</b>	<b>14</b>
<b>Dependant Values: .....</b>	<b>15</b>
<i>The System button</i> .....	16
System Tab .....	16
System Errors .....	16
Windows Event Viewer .....	17
Leak/Compliance .....	17
Piston Initialisation .....	17
Leak & Compliance .....	17
<b>Running the Leak &amp; Compliance procedure .....</b>	<b>18</b>
<b>The Patient Data button.....</b>	<b>21</b>
<i>Data Logging</i> .....	21
<i>Finalise Recording</i> .....	21
<i>Data Storage Locations</i> .....	21
<b>Viewing Recorded Data .....</b>	<b>22</b>
<i>Viewing Recorded data during a recording</i> .....	22
Viewing the instantaneous value using Data View .....	22

# USING THE TAFONIUS SOFTWARE

<i>Viewing Recorded data from the Tafoneius Shell</i> .....	23
<b>Exporting Data</b> .....	<b>23</b>
<b>The Configure button</b> .....	<b>24</b>
<i>Presets</i> .....	24
<i>Initial Settings</i> .....	25
<i>User settings</i> .....	25
<i>Changing Presets</i> .....	26
Save As New .....	26
Save Active .....	26
Rename .....	26
Switch To .....	26
Delete .....	26
Reset Factory Default .....	26
<i>Gas Monitoring</i> .....	27
Gas monitoring tab with the IRMA Gas Analyser: .....	27
Gas monitoring tab with the Artema AION Gas Analyser: .....	28
<i>IRMA Setup</i> .....	29
<i>Pre-Use check</i> .....	31
<i>IRMA Monitor: General Description</i> .....	31
<i>Calibration, self-test and reset</i> .....	31
Oxygen Recalibration .....	31
CO <sub>2</sub> , Nitrous and Agent Calibration (Zero reference calibration) .....	32
<i>Using the ARTEMA AION Gas Analyser</i> .....	33
<i>Display of Measured Values</i> .....	33
<b>Other Control buttons: Alarm Silence &amp; Freeze Traces</b> .....	<b>35</b>
<i>Alarm Silence</i> .....	35
<i>Freeze Traces</i> .....	35
<b>The countdown time for Freeze Traces is determined by the Countdown value set in the Initial settings attribute panel.</b> .....	<b>35</b>
<b>Monitoring Software</b> .....	<b>36</b>
<i>Monitoring options</i> .....	36
<i>Organisation of traces</i> .....	37
Adding and Removing Traces .....	38
<i>Controlling the appearance of traces</i> .....	38
The ECG/IBP Attribute tab .....	40
The CO <sub>2</sub> Attribute tab .....	43
The Airway Pressure (AP) Attribute tab .....	44
The Oxygen Attribute tab .....	45




## USING THE Tafoneius SOFTWARE

The SPO2 Attribute tab.....	46
The Agent Attribute tab.....	47
The N2O Attribute tab.....	48
The Temperature Attribute tab.....	49
The Traces/Trends Attribute tab.....	50
<b>Alarms, Alerts and Critical Error Messages.....</b>	<b>51</b>
<i>Creating an Anaesthetic Record Chart</i> .....	52
Altering the Start Time of the Anaesthetic Record chart.....	52
<i>Viewing the Anaesthetic Record chart</i> .....	52
Data Recording using the Slow Data File.....	53
Changing the interval for Periodic Data in the Slow Data File.....	53
Effect of the Logging Interval on the Slow Data entries.....	53
<b>The Tafoneius Shell.....</b>	<b>56</b>
<i>Tafoneius</i> .....	56
<i>View Data</i> .....	56
<i>Export Data</i> .....	58
<i>Touch Screen Alignment</i> .....	59
<i>Recover Anaesthetic Records</i> .....	59
<i>Company Name</i> .....	59
<i>Shutdown</i> .....	59
<b>Appendix A : SLOW DATA INFORMATION.....</b>	<b>60</b>

# USING THE Tafonius SOFTWARE

## Pre-filling the system

After **Tafonius** has been initialised the piston is left at the bottom of the cylinder ready for filling. Make sure the Y-piece is capped, turn the vaporizer to the required pre-fill concentration and turn on the manual oxygen flow meter to fill the cylinder with gas. The cylinder icon on the left of the screen indicates the volume of gas in litres in the cylinder.



Fill the cylinder to a level equal to the maximum expected patient's tidal volume PLUS 2.5 Litres. e.g. 500kg horse. Average tidal volume predicted is 5.0L so a larger breath of 7.5L can reasonably be expected. Then add 2.5L, which is the buffer volume used by the ventilator. Total pre-fill volume = 10L. When the cylinder graphic shows it is at the required volume, turn off the oxygen flow.

This sets the volume of the “Virtual Bag” for spontaneous breathing. The concept of the Virtual Bag is central to how **Tafonius** operates. This action of the ventilator maintains a minimum system volume thereby reducing anaesthetic and oxygen consumption as well as decreasing the system's time constant.

### The Piston as a “Virtual bag”

When the ventilator is in standby mode and a patient is connected, the system behaves like a ‘perfect’ bag. When the patient breathes out there is no expiratory resistance and when the patient breathes in there is no inspiratory resistance. This ideal performance is accomplished with the Airway Servo System employed by Tafonius. The airway pressure is sensed at the centre of the Y-piece. As soon as a pressure deviation in excess of 0.5cm H<sub>2</sub>O from the ambient pressure is detected then the piston is moved so as to negate this pressure change. If the patient breathes out then the piston withdraws upwards, effectively filling “the bag”. If the patient inspires then the piston moves down, effectively emptying. Because the piston is driven by a motor it is the motor that overcomes the resistances of any piping, valve or soda lime. To the patient it feels like breathing into room air albeit through and ET tube.

### Dictating the size of the “Virtual Bag”

With varying sizes of patients it is beneficial to keep the size of any reservoir bag in the system to a minimum. This then allows rapid responses to changes in anaesthetic or oxygen concentrations. The setting of the Tidal Volume and the concept of fixed buffer volumes control the size of the “Virtual Bag” in **Tafonius**.

tafonius employs two fixed buffer volumes; an upper and a lower. The lower buffer volume is 1.5L and the upper buffer volume is 1.0L. These volumes are constant and are added to the set Tidal Volume to determine the virtual “Bag” size thereby creating upper and lower volume limits.

The upper volume limit = Lower Buffer + Tidal Volume + Upper Buffer.

The lower volume limit = Lower Buffer

Whenever the piston volume exceeds the upper limit, gas is removed from the system. When the piston volume falls below the lower limit then gas is added to the system.

# USING THE Tafonius SOFTWARE

## Starting a Session

Tafonius uses the concept of anaesthetic or procedure **Sessions** for the collection and recording of data. As soon as the Tafonius software starts then information begins to be stored in the Slow Data file and the Session is started. This storage of data continues until a session is Finalised. A session may be Finalised by one of the following:

- Clicking or touching the Finalise Recording button in the Data Logging tab
- Exiting the program through the System and then Exit buttons
- Exiting the program by pressing and holding the Mains Control button
- An abrupt program closure such as power loss, closure of Windows without closing Tafonius (non-controlled shutdown)

In all but the last procedure, the data will be saved **and** associated with all the patient information.

In the event of a non-controlled shutdown there will be NO DATA LOSS, although the patient information will not be associated with the data. In this instance the Files associated with the session will be saved with the word "Patients" as a prefix and the Slow Data and Anaesthetic Record will use only the date as identifiers. For example, a session interrupted by a non-controlled shutdown on January 7<sup>th</sup> 2009 will have a record named "Patients2009Jan07" with Slow Data and Anaesthetic Records named AR2009Jan7-8\_50\_56 and SD2009Jan7-8\_50\_56 respectively. The suffix 8\_50\_56 indicates that the session began at 08:50:56 on the morning of January 7<sup>th</sup>.


For more information on entering Patient and case information, see the section entitled "Entering Initial Patient Information" later in this manual

For more information on the production of an Anaesthetic Record see the section entitled "Creating An Anaesthetic Record Chart" later in this manual

# USING THE TAFONIUS SOFTWARE

## The Finalise Recording procedure

Using any of the first three Finalising methods listed above will cause the following dialog to appear



**Accept Findings**

Do you accept these findings:

Hospital Number:	12345678
Patient Name:	Alfonzo
Owner Name:	Smithers
Patient ID:	192837465
Species:	Equine
Breed/Type:	Arab
Age:	4
Sex:	M
Body Weight:	500Kg
Notes:	Alfonzo notes

At Presentation:

Anaesthetist:	John Irving (House Surgeon)		
Surgeon:			
Procedure:			
Respiration Rate:	0 Breaths/min	Pulse:	0 Beats/min
Temperature:	0 Deg C	CRT:	1sec
PCV:		MM:	Pink
TTP:		Pant:	No
BUN:		Fasted:	No
Anaesthetic Risk Status:	1		

To accept these findings press Accept, to clear the finding press Discard, or to return to the application and change the findings press Change.

This dialog presents the current information associated with the recording session. If these are all correct then click or touch the Acept button. To edit any feature shown, click or touch the Change button. This will take you to back to the Patient Data edit region, where any of the details can be changed.


If all of the details are correct, then clicking or touching the Acept button will cause the following to happen automatically:

- The Session is ended.
- The patient details are used to create the appropriately named Slow Data and Anaesthetic Record files.
- A new Session is started using, either the next set of patient information in the Case list, or the default Patient/Anonymous details.
- Recording for the new Session is started

# USING THE Tafonius SOFTWARE

## Entering initial patient information

Touch the Patient Data button



A new window appears in the information area

A screenshot of the 'Patient Details' window. It has tabs at the top: 'Cases', 'Patient Details' (selected), 'Notes', 'At Presentation', 'Data Logging', and 'Data View'. The form contains fields for: 'Hospital Number:', 'Patient Name:' (with 'Anonymous' entered), 'Owner Name:', 'Patient ID (PID):' (with '111111' entered), 'Species:' (with 'Unknown' entered), 'Breed/Type:', 'Age:' (with '0' in 'y' and '0' in 'm' boxes), 'Sex:', and 'Weight:' (with '0' in 'Kg' box). A 'Close' button is at the bottom right.

Enter the Hospital Number (if used) and then use the tab key to move to the next data entry field. Fields that take numerical data entry will bring up the small on-screen numerical key-pad which can be used instead of a keyboard. By default the patient's name is "Anonymous". Edit this to enter your own patient's name. When the details have been entered touch the "At Presentation" tab.

A screenshot of the 'At Presentation' window. It has tabs at the top: 'Cases', 'Patient Details', 'Notes', 'At Presentation' (selected), 'Data Logging', and 'Data View'. The form contains fields for: 'Anaesthetist:', 'Surgeon:', 'Procedure:', 'Fasted' (radio buttons for Yes/No, with 'No' selected), 'Pant' (radio buttons for Yes/No, with 'No' selected), 'Anesthetic Risk Status' (radio buttons for 1, 2, 3, 4, 5, and E, with '1' selected), 'Temp' (radio buttons for Fahrenheit/Celcius, with 'Celcius' selected), 'RR:', 'PCV:', 'TP:', 'BUN:', 'Pulse:', 'CRT:' (with '1sec' selected), and 'MM:' (with 'Pink' selected). A 'Close' button is at the bottom right.

Fill in the procedure details along with the rest of the presenting clinical information.

## USING THE TAFONIUS SOFTWARE

Entries can be left blank. The Anaesthetist, Surgeon and Procedure fields are used when an Anaesthetic Record chart is produced. The Anaesthetist, Surgeon and Procedure entries are all taken from a Drop-Down list. To view the list click or touch the down arrow in the blue box next to the field you wish to enter details for. If the item you require appears on the list click or touch it to select it. This provides a quick way to enter commonly used information. If the item does not appear on the list then it may be added by touching or clicking on the ellipse button



Clicking or touching the Ellipse button opens an editor.  
Here the Anaesthetist Drop-Down list has been opened for editing:

Anaesthetists	
Dr Andrews	
Lorraine Gosch	
Dr Simpson	
John Irving (House Surgeon)	

Items in the list may be ordered using the Up & Down buttons. The Add, Edit and Delete buttons allow standard editing of the entries in the list. When editing is complete, click or touch the Close button to return to the “At Presentation” tab.

### Notes Tab:

Any notes that you want to enter relating to the case can be entered as free text in the Notes tab. These notes will appear in the Anaesthetic Record chart in the area under the Header of general information.

## USING THE Tafonius SOFTWARE

### Cases Tab:

Details relating to cases can be pre-entered as a list of patients. Associated case data can be entered and then recalled when required.

Patient Name	Patient ID (PID)
Anonymous	0000000
Alfonzo	192837465
→ Diamond Sunset	198725463

Select

New

Delete

By default, when the software is started details appear for a patient called Anonymous with a patient ID of 000000.

The green arrow next to the patient name shows that these patient details are in use. Patient details cannot be deleted if they are in use: you must first select another patient. Touching or clicking the New button opens the Patient details tab ready for new patient information. To change patients, click or touch the patient name or ID and then click or touch the Select button. The selected patient information is then used for the current recording session. When the monitoring/recording session is finalised, the Patient details in use at that time will be the ones used to create the recording files. During the Finalising procedure a dialog box will appear to check that the Patient details to be saved are correct.

The patient name and ID is always shown at the bottom of the Tafonius screen. Click or touch either of these two screen areas to open the Patient attribute dialog.

## USING THE Tafonius SOFTWARE

### Data Logging Tab:

In this tab the options for Slow Data Recording and Anaesthetic Record chart heart rate can be set. Because Tafonius has 3 means of monitoring heart rate (ECG, IBP & Pulse-Ox) it is necessary to select one to be used for the Anaesthetic Record chart. Unless you have specific reasons for changing the Slow Data Logging time, leave this set at 5 seconds.

The screenshot shows the 'Data Logging' tab selected in the software interface. The tab bar at the top includes 'Cases', 'Patient Details', 'Notes', 'At Presentation', 'Data Logging' (highlighted), and 'Data View'. The main content area contains the following elements:

- Slow Data Logging:** A text label followed by a numeric input field containing '5' and the unit 'sec'.
- Time:** A label in blue text.
- Start:** 14:27:55
- Elapsed:** 00:04:14
- Record Spontaneous Breaths:** An unchecked checkbox.
- Record Ventilated Breaths:** A checked checkbox.
- For Anaesthesia Record Log Heart Rate From:** A section with three radio button options: 'ECG' (selected), 'IBP', and 'Pulse-ox'.
- Anaesthetic Record Start:** A button.
- Finalise Recording:** A button.
- Close:** A button.

When all the patient information has been entered, click or touch the “Close” button.




# USING THE TAFONIUS SOFTWARE

## The Ventilator controls

The ventilator controls consist of Tidal Volume, Respiratory Rate, Inspiratory Time, Maximum Working Pressure Limit and a CPAP (or PEEP) setting. These are the only controls that can be set for ventilation.

The main controls (Tidal Volume, Respiratory Rate & Inspiratory Time) have two sets of numbers associated with them. These numbers change depending on whether the system is in Standby (spontaneous breathing) mode or Ventilate mode. In Standby mode, the large blue numbers are the measured values derived from the patient's respiratory efforts. The Airway Servo System is working as a respiratory monitor. The smaller black numbers are the ventilator settings that will be used in Ventilate mode. If no breathing is detected after 60 seconds the blue values are all set to zero.



When you touch a button it changes colour and becomes ready to accept a new value. A numeric keypad appears in the lower right hand corner of the screen where you can enter a new value. Or the value may be changed by means of the UP/DOWN arrows. When the new value has been entered do one of the following:

- Touch the body of the button again,
- Touch the "Enter" button on the numeric keypad or
- Touch the body of the next button whose value you want to change.

## USING THE Tafonius SOFTWARE

### Tidal Volume:

In Standby mode: The set Tidal Volume is indicated by the smaller black number and is the volume that will be delivered per breath when the switch is made to Ventilate mode. It also sets the volume of the Virtual Bag used by the patient. The blue value indicates the breath-by-breath tidal volume of the spontaneously breathing patient

In Ventilate Mode: Sets the actual volume delivered to the patient. **Tafonius** uses automatic compliance compensation so the volume received by the patient is the volume entered. The volume displaced by the cylinder will be more than the set value to compensate for the volume lost to the compliance of the system. The patient actually receives the TV you set regardless of the length of the breathing systems tubes. This is very helpful when 20 or 30 ft tubes are connected to reach a patient in an MRI. Tidal Volume units are Litres

### Respiratory Rate:

In Standby mode: The blue figures show the measured respiratory rate, recalculated every 10 seconds on a rolling 1 minute basis. This gives a true reflection of respiratory rate. The smaller black figure shows the respiratory rate that will be used in Ventilate mode.

In Ventilate Mode: The large black number shows the respiratory rate in use.

Respiratory Rate units are Breaths per Minute

### Inspiratory Time:

The Inspiratory Time button controls the total length of the Inspiratory phase. There is also a smaller sub-button located within the Inspiratory Time button that allows a pause to be included in the Inspiratory time. This inspiratory pause is limited to 40% of the I time. It can be displayed in either seconds or percent of I time, by repeated touches of the sub-button.

In Standby Mode: The blue figures show the measured inspiratory time, recalculated after every breath i.e. on a breath to breath basis. The smaller black figure shows the total Inspiratory Time that will be used in Ventilate mode.

In Ventilate Mode: The large black figure indicates the total time for inspiration. This may not be the same as the Inspiratory flow time, depending on the setting for Inspiratory Pause.

### Inspiratory Pause:

In Standby Mode: Has no action

In Ventilate Mode: Controls the portion of the Inspiratory phase that is held as a pause. The maximum pause is limited to 40% of the total Inspiratory time. Increasing the Pause time will increase the time the inspiratory volume is held at the end of the Inspiration stroke. The total Inspiratory Time is unaffected so Inspiratory Flow rates will increase as more Pause time is added.

Touch the Pause button to activate it and open the numeric keypad. Touch it again to toggle between setting the pause as a % of Inspiratory time, or as a finite time. Enter the actual time or percentage required. Press "Enter" on the numeric keypad when done.

**In general, as with the ventilator settings above, any button shaded green on the display is a value you can set. These include the following.**

## USING THE Tefonius SOFTWARE

### Maximum Working Pressure Limit (MWPL)

The Maximum Working Pressure Limit control sets a level, above which ventilation is prevented. This is particularly useful when dealing with neonates or during open-chest surgery when lung expansion is not restricted by a closed thorax. The Maximum Working Pressure Limit is an absolute value and is unaffected by any settings of PEEP or CPAP. If the MWPL value is reached during an Inspiratory phase then the ventilator immediately sounds an alarm and changes to the expiratory phase, where once again the cylinder behaves as a Virtual Bag. The MWPL setting can be useful for patients that are “fighting” the ventilator to prevent excessive pressures developing.

### CPAP/PEEP

In Standby Mode: CPAP or Continuous Positive Airway Pressure can be applied during spontaneous breathing. . With CPAP patients inspire and expire normally with no added effort, although the end expiratory pressure is elevated to the CPAP setting. This aids in maintaining open alveolae. The airway pressure is held constant during all phases of breathing.

In Ventilator Mode: PEEP or Positive End Expiratory Pressure applies during IPPV. Ventilating pressures rise during the Inspiratory phase as normal, but are allowed only to fall as far as the PEEP setting during expiration

### Validating entries

All entries are validated once complete. This means that any invalid entries that would exceed minimum or maximum settings that **Tefonius** is capable of delivering will not be accepted. When this occurs a warning sound is issued and a description of the input error is shown at the bottom of the screen as a black on yellow text message. In most circumstances the limiting factor preventing the input is also shown. Here, for example a value of 1.0 second has been entered as a Pause Time during an Inspiration Time of 2.0 seconds. This represents 50% of the Inspiratory Time, thus exceeding the maximum of 40% and so is disallowed. The maximum allowable entry of 0.8 seconds is used.

Increasing the It to 2.5 seconds would allow a 0.5 second pause. The Pause Time box is outlined in red to indicate the source of the error.

Tidal Volume 5 <b>0</b> (Litres)	Respiratory Rate 5 <b>0</b> (Breaths per Minute)	Inspiratory Time 2 <b>0</b> (seconds)
Minute Volume <b>0</b> 25 (litres)	Inspiratory Flow <b>0</b> 250 (litres/minute)	I:E Ratio 1:5 <b>1:0</b>
Expiratory Time <b>0</b> 10 (seconds)		

The Inspiratory Pause Time setting 1s was out of range. Valid range is 0-0.8s. Change limited by IT.

## USING THE Tafonius SOFTWARE

### **Dependant Values:**

The four buttons below the ventilator controls show the derived values for Minute Volume, Inspiratory Flow, I:E Ratio and Expiratory Time. They are dependent upon the settings of TV, RR and Itime


In Standby Mode: The large blue figures show the values calculated from the measurements made during normal inspiration and expiration. Respiratory rate is calculated on a rolling minute average that is updated every 10 seconds. The measured minute volume is an integral of the ventilated volume over the last minute. Again this is updated every 10 seconds. If no respiration is detected during a period of 60 seconds then these blue values are all reset to zero. The smaller black figures show the calculated values based on the values set for the three independent ventilator settings (TV, RR & IT).

In Ventilate Mode: The large black figures show the calculated values based on the values set for the 3 independent controls (TV, RR & IT).

# USING THE Tafonius SOFTWARE


## Control buttons: Controlling the Software

There are three buttons that control the look, feel and behaviour of the **Tafonius** software. These three buttons – System, Patient Data & Configure are located at the bottom right hand side of the screen.



### The System button

The System button opens a two-tab screen that contains a System and a Configuration tab. The Configuration tab can be entered directly from the “Configure” button on the lower right hand side of the screen.



### System Tab

The System tab contains buttons for System Errors, Leak/Compliance test settings, Language and User Manual options. Currently the Language is set to English only. Other language options will become available in due course. The User Manual button also presents options which at the time of writing are not yet available.

### System Errors

Any error or fault condition that occurs in the **Tafonius** software is reported in this System Errors list. An error is indicated by the flashing of the border of the System button in red and an audible alarm. The purpose of the list is to indicate the source of an unexpected alarm or fault condition and to assist with trouble-shooting. Some errors are transient and will self-clear, in which case no action is

## USING THE Tafonius SOFTWARE

necessary. Some errors will not self-clear or will remain for example after a device has been removed. In these instances, select the error and then click or touch on the “Resolve” button. If the error is no longer present then it is resolved and disappears from the list.

### Windows Event Viewer

This is reserved for use as a trouble-shooting tool. A field service engineer or someone from technical support may ask you to use this button to view any error sequences that have arisen during the use of Tafonius. Do not use this button during normal use as there may be an effect on the performance and behaviour of the Tafonius software.


### Leak/Compliance

A leak test and Compliance measurement is offered whenever the **Tafonius** software is launched. However sometimes it is convenient or necessary to run either utility without having to close the program and start again. Clicking on the Leak/Compliance button starts a sequence of events that proceeds in the following order:

(If the piston is not initialised then Piston Initialisation must be done before either leak testing or compliance measurement.)

### Piston Initialisation

This involves making sure that there is nothing attached to the Y-piece before the piston is sent to the bottom of the cylinder to identify and reset its zero position. If the piston has not been initialised since turn-on then the cylinder representation on the screen will have a hatched background. Once the piston has been initialised this changes to a solid green background. Leak testing and compliance values can then be measured.



Uninitialised piston background is hatched. Image also has a red flashing border

Initialised piston background is green and shows piston volume in litres




### Leak & Compliance

This button begins the process of first leak testing and then measuring the compliance of the system. For leak testing the piston is taken to the mid-point of the cylinder and then the piston produces a pressure of 20cm H<sub>2</sub>O. During this pressure phase the piston movement and hence leak is measured. The measured leak should be below 1500mls/min and in normal circumstances is less than 500 mls/minute. Values above or around 1500 mls/minute represent a serious leak and should be investigated. However, due to the way **Tafonius** maintains pressures and controls piston movement, even leaks of 1500mls/minute are unlikely to affect ventilator performance. They do however represent a loss of gas and agent.


# USING THE TAFONIUS SOFTWARE

## Running the Leak & Compliance procedure


Press the Leak/Compliance button on the System tab to begin the procedure




A dialog appears announcing that it is establishing a connection with the Auxiliary controller. This can take 5-10 seconds. After that a second dialog appears asking that any patient be removed from the circuit. When the breathing circuit is ready, press the Leak & Compliance button



At this point a stopper must be applied to the Y-piece to seal the system. Then press the Leak & Compliance button.




The piston is moved to the mid-point and then stops briefly before establishing a system pressure of 20cm H2O.




## USING THE Tafarius SOFTWARE

If for any reason there is an excessive leak the following message appears:



In this instance seal the leak and press the Retry button.


If the leak test is successful the leak value is displayed and the system then proceeds to perform the two compliance measurements.



A compliance measurement is made at the top cylinder position and at the 5.0L cylinder position and uses the leak value obtained in the first test. At the end of the test the leak and both compliance values are shown. Note that the values obtained will vary markedly depending on the nature of the breathing circuit connected. The two compliance values are used to obtain the equation for the straight-line graph of compliance for your machine and the breathing system in use at that time. This equation is then used to compute the compliance compensation required at any piston position and pressure.



## USING THE Tefonius SOFTWARE



Press "Finish" to leave the Leak & Compliance test procedure.

# USING THE **Tafonius** SOFTWARE

## **The Patient Data button**

Entering patient data has been described under the section “Entering Initial Patient Information”. This covers the “Patient Details” and “At Presentation” tabs found once the “Patient Data” button is pressed. Also found here are the “Data Logging” and “Data View” tabs.

### **Data Logging**

All physiological data and many automatic events are recorded in a file called the Slow Data (SD) file. How often data is entered in this SD file is dictated by the setting in the Data Logging tab. The default value is 5 seconds. A range between 1 second and several minutes can be entered.

Some natural events such as inspiration and expiration can be entered into the SD file. This however can fill the file with information that may not necessarily be required. Spontaneous breaths or ventilated breaths can be selected or deselected for entry in the SD file by checking or un-checking the respective tick-boxes.

The time that a particular recording session began and its duration are also shown in this tab.

### **Finalise Recording**

As soon as the **Tafonius** software starts, data is recorded in the SD file. When a patient session is ended the recording should be finalised. To do this push on the “Finalise Recording” button. A dialog box may appear asking to save the current Preset (see the section on User Presets) and then a further dialog box will appear stating that the data will be committed and that no further changes can be made. Click or touch OK to commit the data and begin a new session. When a session is Finalised, those patient details are removed from the Patient List in the Cases tab.

### **Data Storage Locations**

Where the data is stored will depend whether any Patient Information has been entered.

If no Patient data has been entered then the data is stored in a folder that is named using the word “Patients”, the calendar month and the day. Several patients or sessions may be held in this folder if no patient information is entered during any of the day’s operations.

e.g. An operation on September 12<sup>th</sup> 2008 where no patient information was entered will be stored in a folder called “Patients2008Sep12”

If patient data has been entered then the folder name will be derived from the Patient’s name and the patients ID.

e.g. Wild Billy , ID 1786394 will have a session folder created called WildBilly\_1786394. If there is more than one session recorded for Wild Billy on that day then all sessions will be stored in this folder.

Recorded data can be easily viewed and or exported using the utilities supplied with the **Tafonius** software. See the sections “Viewing Recorded Data” and “Exporting Data”


# USING THE TAFONIUS SOFTWARE

## Viewing Recorded Data


Recorded data can be viewed by two methods: During a recording session in the ventilator control software or after a session in a spreadsheet viewer.

### Viewing Recorded data during a recording

To view data on-the-fly, click or touch on “Patient Data” and then the “Data View” tab.




A trend viewer pops up showing the trend of the selected parameter from the beginning of the recording session to the current time. This data is constantly updated at a rate set by the Slow Data logging period. Use the drop-down list to select a different parameter.



### Viewing the instantaneous value using Data Drop Down List

Click or touch a point in the graph area of the Data View tab. A line appears at the point of touch. At the bottom of the line on the x-axis of the graph is the time at which that point of data was recorded. Along the right hand side of the line are the instantaneous value(s) of the associated data.



## USING THE **Tafonius** SOFTWARE

### **Viewing Recorded data from the **Tafonius** Shell**

Click or touch on the “View Data” button in the **Tafonius** Shell screen. This opens a Record Picker box to choose the Data Record to view. Data is presented in date order with the most recent at the top. Pick a record to view and then click/touch the “Select” button. A file or list of files will be shown that begin with SD. These are the Slow Data files recorded during a session. Touch or click or touch on the one required and then press the “Select” button. A spreadsheet viewer opens to display the recorded data. Touching a cell within a column representing measured physiological data will display that data automatically as a trend graph.

To create graphs for other parameters that are not physiological such as the set Tidal Volume or Compliance Compensated Tidal Volume it is necessary to export the data to a spreadsheet viewer of your choice. This is easily done from the **Tafonius** Shell.

### **Exporting Data**

From the **Tafonius** Shell, click or touch on the “Export Data” button. An “Export Data” dialog box opens listing the recorded data on the **Tafonius** computer. Data is presented in date order with the most recent at the top. Choose the Record folder required and then click or touch on the “Next” button. Then choose the destination for your exported data. Typically this will be a USB flash drive inserted in the hub at the front of the screen. When the destination has been set, click or touch on the “Export” button to copy the data across to the selected device. Data is exported as a Comma Separated Value (CSV) file. The Anaesthetic Record chart in pdf format will be exported at the same time.

# USING THE Tafonius SOFTWARE

## The Configure button


The “Configure” button is used to access the following: Presets, Initial Settings, Gas Monitoring and SPO2 setup. The “Security” button in the Configure window is not active at this time.

### Presets

Presets are a means of individually saving **Tafonius** software layouts and initial settings that suit a particular purpose. For example you may have a preset called “Foals” or a preset called “Colics” or even “Senior Anaesthetist”.

Each Preset saves the way the screen is configured including how the traces are arranged as well as all the initial values used when the program starts, as well as all the alarm settings and whether they are enabled or not.

There are two ways in which Presets are used to load values at the start of a procedure. There are the Initial Settings, which describe the preferred initial values for such items as Tidal Volume, Respiratory rate etc and there are the user settings which remember the preferred settings for e.g. ECG speed and trace position. The values for Initial Settings can be changed by clicking on the Configure button at the bottom right hand side of the screen. A new smaller window area appears below the traces area.



Choose the preset you want to change by clicking on it in the Presets list. If it is already in use it will have the word Active after the preset name. If not, then click or touch on the required Preset and then choose the “Switch To” button. The selected Preset will be loaded. Then click or touch on the Initial Settings button on the left hand side.

# USING THE Tafonius SOFTWARE

## Initial Settings

Using the keyboard or the on-screen keypad enter the **initial** values you would like in each of the named fields.

The fields have the following meanings:

**TV:** Tidal Volume in litres.

**RR:** Respiratory rate in breaths per minute

**It:** Total Inspiratory time

**MWPL:** Maximum Working Pressure Limit

**IP:** Inspiratory Pause in seconds

Show IP as % ☐ : Shows the Inspiratory pause as a percentage of the **It**.

**CPAP:** Continuous Positive Airway Pressure in cm H<sub>2</sub>O

**FiO2:** Desired initial oxygen fraction. (not available in software version 1.0)

**APNOEA:** Length of time allowed between breaths before Apnoea alarm sounds

**Countdown:** Length of countdown timer before pop-up menus close automatically, in seconds

**Countdown** ☐: Enable the countdown timer for closing pop up menus

**Mouse** ☐: Enable the mouse cursor on the screen

When done, click or touch on the Preset tab again and then choose "Save Active". The next time this Preset is loaded these new initial values will be used.

## User settings

During the course of using the **Tafonius** program you may make changes to the screen appearance or to the appearance of a certain monitor trace. When changes such as these are made, that affect the personal way that the screen is viewed the Preset will be marked as altered and an asterisk will appear next to the name in the title bar.



At any time or at the end of a procedure this unique set up of Tafonius can be saved by again choosing to save the Preset. If this is not done explicitly during use then you will be prompted to save the Preset when finalising a patient session. You can choose to change the active preset by touching SAVE or you may save the current setup as a new setup by touching SAVE AS NEW and renaming it as desired.


Note that you will not be prompted to save a changed Preset when you shut down the Program using the main ON/OFF switch.

# USING THE Tefonius SOFTWARE

## Changing Presets

When the **Tefonius** software starts it uses, by default, the last Preset in use. To change a Preset, click or touch on "Configure". The Preset list is displayed. Along the bottom of the Preset Window are 5 buttons, used to control how Presets are created, named, selected and deleted plus two other buttons to Reset Factory Defaults or Close the Configuration dialog.

The Preset currently in use is indicated by the word Active in brackets after the Preset name. This name is also displayed at the top of the screen in the title bar.



### Save As New

This button creates a copy of the Active Preset and names it "Custom Preset". Use the Rename button to change the name to one of your choice. If you have made changes to the stored active preset the new stored preset will include those changes.

### Save Active

This button saves any changes made to the Active Preset

### Rename

This button allows a Preset to be renamed. Type the new name in the edit box that appears and click or touch "OK".

### Switch To

Use this button to change to a different Preset from the Active one. To use this button, click or touch on the Preset you wish to use and then click or touch on the "Switch To" button.

### Delete

This button removes the highlighted Preset. Note that it is not possible to delete an active preset.

### Reset Factory Default

This button will apply all the default factory settings to the Active Preset. Note that the Preset must be Active for the changes to be made to it. A confirmation dialog box appears stating which Preset will be changed. If this is the Preset that you wish to Reset, then click or touch "OK".


# USING THE TAFONIUS SOFTWARE

## Gas Monitoring

This tab displays all of the status and diagnostic information for the current Gas Monitor in use. At the current time a mainstream device called IRMA or a sidestream device called Artema AION is used for CO<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>O and Agent monitoring. Your machine will be fitted with one or other of these units.

### Gas monitoring tab with the IRMA Gas Analyser:

Clicking on the Gas Monitoring tab reveals a Window with 5 further tab options: Status, Sensor, Errors, Config & Data



### Status tab

Used for fault diagnostic information only

### Sensor tab

Used for fault diagnostic information only

### Errors tab

Used for fault diagnostic information only

### Config tab

Used for fault diagnostic information only

### Data tab

Location of the Zero option for the IRMA monitor. See the section on "Use of the IRMA Monitor" for full information on when and how to use the "Zero" option found in this tab.

In most situations it is only necessary for users to use the Data tab. Information from the other tabs may be requested by a service engineer or a technical support person during fault diagnosis.



## USING THE TAFONIUS SOFTWARE

Gas monitoring tab with the Artema AION Gas Analyser:

The screenshot shows the Tafonius software interface. On the left is a vertical sidebar with buttons: 'Presets', 'Initial Settings', 'Gas Monitoring' (highlighted in yellow), 'SpO2', and 'Security'. The main area has a top bar with 'Configuration' and 'System' tabs. Below this is a sub-bar with 'Modes/Status', 'Status', 'Alarms', 'Cmds', 'Demo', and 'Calib' tabs. The 'Modes/Status' tab is active, displaying a 'Reset' button, 'Op Mode: ISOAccuracy', 'Comm Mode: Normal', 'Atm (hPa): 102.1', 'Water Trap: Adult', 'Pump: Regulating', 'Flow: 202', and a 'Message: Command OK' field. A 'Close' button is in the bottom right corner.

### Modes/Status

Used for fault diagnostic information only

### Status

Used for fault diagnostic information only

### Alarms

Used for fault diagnostic information only

### Cmds

Used for fault diagnostic information only

### Demo

Used for demonstration purposes only.

### Calib

Used for fault diagnostic information only

In most situations it is only necessary for users to use the Cmds tab for performing a zero calibration. Information from the other tabs may be requested by a service engineer or a technical support person during fault diagnosis.

## USING THE Tafoneius SOFTWARE

### Using the PHASEIN IRMA Monitor

#### Important Note:

*The IRMA mainstream multi-gas probe is intended to be connected to a patient breathing circuit for monitoring of inspired/expired gases of patients in intensive care, anaesthesia and emergency care. It is not intended to be used as the only means of monitoring a patient. It shall always be used in combination with other vital signs monitoring devices and/or professional human judgements of patient condition.*

*The IRMA multi-gas analyser is intended to be used by trained and authorised veterinary professionals only. It is not intended to be used in outdoor transport applications such as in cars or in aircrafts.*



#### IRMA Setup

**Plug the IRMA device into the IRMA socket located on the bulkhead behind and below the utility drawer, on the Right Hand Side of the machine.**

**Snap the IRMA sensor head on top of the IRMA airway adaptor. It will click into place when properly seated**



**A green LED indicates that the IRMA sensor is ready to use**



## USING THE Tefonius SOFTWARE

Connect the 15mm male IRMA/Airway adaptor connector to the breathing circuit. When used with a foal 22mm circuit connect the IRMA device as shown in the next 2 images.



Connect the 15mm female IRMA/Airway adaptor connector to the patient's endotracheal tube



When used with a 50mm mainstream adaptor the IRMA adaptor is integrated into the Large airway adaptor as shown:



## USING THE Tafonius SOFTWARE

**Always position the IRMA sensor with the O2 cell (or top) pointing upwards**



### Pre-Use check

Before connecting the IRMA adaptor to the breathing circuit, verify the O2 calibration by checking that the O2 reading on the monitor is correct (21%). See the following section “**Oxygen Recalibration**” on how to perform room air calibration.

### IRMA Monitor: General Description

The IRMA breath monitor unit is a small cube-shaped unit housing a replaceable Airway Adaptor and Oxygen sensor. The main housing for the unit has an LED in one corner which indicates the status of the sensor. These states are described in detail in the section "IRMA Status & Error Messages". The IRMA unit connects to the **Tafonius** monitoring system via a 9-Way Female D Type connector. When the unit is connected and the unit powered ON, the IRMA device runs through a self-test sequence. This is indicated by the LED on the main housing changing alternately through GREEN, RED, BLUE a number of times. Once the initialisation sequence is complete the LED stays GREEN. If the IRMA unit has the capability of detecting anaesthetic agents then the LED will turn BLUE in the presence of these agents.

IRMA units may be connected/disconnected while powered without any problems.

### Calibration, self-test and reset

The IRMA unit is self-calibrating and does not require scheduled calibration on a monthly or even annual basis. Calibration is achieved in a number of ways but must be done while the unit is plugged in and powered ON.

**ALL CALIBRATION MUST BE DONE WITH THE UNIT IN FREE AIR AND NOT WHILST CONNECTED IN A PATIENT CIRCUIT**

### Oxygen Recalibration

For calibration of the Oxygen sensor, removal of the Airway adaptor forces an automatic recalibration to 21%. For this reason it is important that the removal and refitting of the airway adaptor is done in conditions where the airway concentration is at 21% such as in a well ventilated room.

1. Remove the Airway Adaptor by pulling the white plastic connector away from the main housing.
2. Wait until the LED flashes RED and the alarm sounds on the monitor.
3. Replace the Airway Adaptor.
4. Check that the LED turns green
5. Check that the O2 reading on the Monitor is 21%

## USING THE TAFONIUS SOFTWARE

### CO<sub>2</sub>, Nitrous and Agent Calibration (Zero reference calibration)

These can all be reset from the Windows software via the Configure/Gas Monitoring/Data tab route. Only perform this function when asked to do so by the Windows program or when indicated by the following conditions:


- 1) A repeated "Gas Concentration out of Range" error message that appears when measuring room air and is not cleared by disconnecting and re-connecting the IRMA device.
- 2) Incorrectly reported high Inspiratory CO<sub>2</sub> concentrations

**ALL CALIBRATION MUST BE DONE WITH THE UNIT IN FREE AIR AND NOT WHILST CONNECTED IN A PATIENT CIRCUIT**

**THE IRMA UNIT MUST BE ALLOWED TO WARM UP AND BE RUNNING FOR AT LEAST 15 MINUTES BEFORE A ZERO REFERENCE CALIBRATION IS PERFORMED**

To RESET the IRMA device, click or touch the "Zero" button in the Data tab of the Gas Monitoring Window found in the Configuration Setup.

A warning dialog box will appear with the following message:



Make sure the conditions for re-calibration are correct before continuing, then click or touch Yes. If you are unsure about the conditions required for calibration, click or touch "No" and read the above paragraphs carefully before running the "Zero" option again.

## USING THE TAFONIUS SOFTWARE

### Using the ARTEMA AION Gas Analyser

The Artema AION Gas Analyser is a sidestream gas analyser, drawing the sampled gas from the Y-piece. Other than the physical differences and the inherent latency between breathing and waveform display, the reporting and display of all gas parameters is the same as with the IRMA Gas Analyser.



Sidestream sampling port (capped)

The sidestream sampling line connects to the sampling catheter that enters the Y-piece as shown. The sampling catheter extends down into the ET tube for optimum sampling.

### Display of Measured Values

THE Gas Analyser (IRMA or AION) will continuously display and update the numerical and trace information for the patient. The trace region shows the capnogram as either a line or filled trace. The main numerical data region associated with the trace shows the inspired and expired values of CO<sub>2</sub>.



The larger figure represents the expired value and the smaller figure shows the inspired value. The third value to the right is the respiratory rate in breaths per minute. If a period of apnoea occurs then the alarm will sound and the respiratory rate will be replaced by a red box, which displays the number of seconds elapsed since the last breath.



## USING THE Tafoneius SOFTWARE

Clicking or touching this red box will silence the Apnoea alarm. The box will remain red and continue to count the seconds since the last recorded breath. A fresh breath will clear the Apnoea state and return to displaying Respiratory Rate (RR).

During periods of Apnoea, the inspired values are constantly updated indicating the gas concentrations in the breathing circuit.

DRAFT

# USING THE TAFONIUS SOFTWARE

## Other Control buttons: Alarm Silence & Freeze Traces

### Alarm Silence

Press the button once for a 5-second silence of all alarms. After 5 seconds the alarms, if active, will again sound.




Press and hold the Alarm Silence button until it shows the legend "Locked".



The audio portion of all current alarms are then suppressed. Any new alarm condition will override the locked out condition and again sound the alarm.

### Freeze Traces

Use this button to freeze all trace information on the screen. All other monitoring activities continue as normal. The button will change colour to indicate that it is active.



Also a countdown timer will begin showing on the button. At the end of the countdown period the freeze option will be lifted. Alternatively simply clicking or touching the Freeze button once more will release the freeze option.

The countdown time for Freeze Traces is determined by the Countdown value set in the Initial settings attribute panel.



# USING THE TAFONIUS SOFTWARE

## Monitoring Software

### Monitoring options

Tafonius has the options of monitoring the following parameters via its in-built software:

Single Channel ECG

Single Channel Invasive (Direct) Blood Pressure

Breath Analysis:

Inspired/Expired CO<sub>2</sub>

Inspired/Expired O<sub>2</sub>

Inspired/Expired N<sub>2</sub>O

Inspired/Expired Anaesthetic Agent:

Halothane

Isoflurane

Sevoflurane

Desflurane

Enflurane

Single Channel Pulse-Oximeter


Dual Channel Temperature with delta

Airway Pressure

# USING THE Tafonius SOFTWARE

## Organisation of traces

The monitoring screen consists of 3 main regions – a trace region where all time-traces are displayed, a numerical region where all values are displayed and a status region where status information about **Tafonius** is displayed.




The trace region is centralised on the screen and in the picture above contains the ECG, CO2, Airway Pressure and Pulse-Oximetry trace. This area can be organised as required by placing any trace or combination of traces here. The trace region itself is divided into three horizontal areas each of which is one trace high. It is possible to place a trace in any of these 3 horizontal areas. Additionally it is possible to have multiple traces in any one of these horizontal areas.

## USING THE Tefonius SOFTWARE

### Adding and Removing Traces

To change the trace screen viewing options touch any trace. A new editing box appears at the bottom of the screen




The text in the yellow strip indicates the trace that has been touched (selected). Here the user has two options: Touching inside the Trace/Trend dialog box will remove the trace from the screen and place it in the reserve area. Alternatively touching the screen in another area of the trace region will place that trace alongside and to the left of any other traces in that horizontal area.

Any monitor that has been placed in the reserve area can be put back on the screen in any of the three horizontal areas at any time. To do this touch any "number box" e.g. the ECG heart rate box. The set of tabbed attribute panels appear at the bottom of the screen. Select the Traces/Trends tab and then select and drop the desired monitor trace back on the screen. "Pick and put" the desired trace where you want it.

### Controlling the appearance of traces

All traces have various attributes associated with them that determine for example their height, speed, shift and alarm options amongst others.

To access these options for any monitor touch the number panel associated with the monitor. For example to access the options for the ECG monitor, touch the left hand side of the ECG/IBP panel at the top right of the screen.



## USING THE TAFONIUS SOFTWARE


The monitoring attribute panels appear at the bottom of the screen opened at the ECG tab. Note that the ECG attribute tab shares its place with the IBP attribute tab as they are both located on this common ECG/IBP tab. When the left hand side of the ECG/IBP panel is touched, the ECG attributes are shown. If the right hand side of the ECG panel is touched, which has the IBP values, then the attribute tab will be opened with the IBP attributes showing. To change between ECG and IBP attributes touch the appropriate radio button at the left hand side of the ECG/IBP tab

All alarm settings and trace settings are particular to any given Preset and so the Preset should be saved after making adjustments to any of the controls in the monitoring attribute panels.

## USING THE TAFONIUS SOFTWARE

### The ECG/IBP Attribute tab

There are two active sections to this tab: Traces and Alarms. In the Traces section all the controls relate to how the trace can be viewed. Values can be set for either ECG or IBP options by selecting the appropriate option in the Select box.



### Trace options: ECG

**Sensitivity:** controls the height of the ECG trace seen on the screen


**Sweep:** Controls the screen sweep speed in rates of mm/second.

**Filter:** Allows the setting of a digital filter to remove unwanted signal noise. Use a notch filter for mains interference. Use a Low-pass filter for tremors due to movement e.g. shivering. Move the slider in the Filter box to control the degree of filtering.

**Shift:** Shifts the base-line of the ECG up and down the screen. Use the shift option to centralise the ECG for optimum viewing.

**Heart Rate:** Select the box next to "Beep on ECG" to produce an audible beep with every detected ECG complex. Since **Tafonius** only has one audio output device, only one monitor can be made to beep at any time. Selecting ECG for example, will mean that any previous selection is overridden.

To change between ECG and IBP attributes touch the appropriate radio button at the left hand side of the ECG/IBP tab



### Synchronising the ECG and IBP waveforms

ECG and IBP waveforms can be synchronised in time as long as they have the same trace width on the screen. Normally they would be placed one above the other so that a visual association between the

## USING THE Tafonius SOFTWARE

two traces can be made. If the two traces are not synchronised then select the “Synchronise” option in the ECG/IBP tab. If the two traces do not have the same width they cannot be synchronised and the option will be greyed.

### Trace options: IBP

The screenshot shows the 'ECG/IBP' tab selected. The 'Trace' button is highlighted. The 'Alarms' button is below it. The 'Select' section has 'IBP' selected. The 'Synchronise' checkbox is checked. The 'Sensitivity' slider is set to 30 mmHg/div. The 'Sweep' section has '12.5 mm/sec' selected. The 'Shift' slider is at 0. The 'Heart Rate' section has 'Beep On IBP' unchecked. The 'Zero IBP' button is visible. The 'Close' button is in the bottom right corner.

### Sensitivity

Moving this slider will change the full scale height of the IBP trace in accordance with the legend on the right hand side of this scale. Note that there is an automatic levelling feature in the IBP trace that means that the bottom of the IBP trace is adjusted after every sweep to lie on the base line. This will always keep the IBP trace in view even if the blood pressure values are slowly changing.

### Sweep

Selecting one of the radio buttons will set the sweep speed for the IBP trace. Note that this option may be further controlled or changed depending on the state of the Synchronise check box and, or the presence of ECG and IBP traces on the Tafonius screen. If the ECG and IBP signals are not synchronised then the sweep speed of ECG and IBP traces can be different. If the traces are synchronised then the IBP sweep speed will be set to and follow the ECG sweep speed.

### Shift

The Shift slider alters the position of the baseline of the IBP trace.

### Heart Rate

Select the box next to “Beep on IBP” to produce an audible beep with every detected pulse waveform. Since **Tafonius** only has one audio output device, only one monitor can be made to beep at any time. Selecting Beep On IBP for example, will mean that any previous beep selection is overridden.

### Zero IBP:

When a direct arterial sampling line has been placed and the sensor positioned at the level of the patient’s heart (but with connection to ambient pressure not arterial pressure) the “Zero IBP” button is used to set the current output to zero.

## USING THE TAFONIUS SOFTWARE

### Alarm options: ECG & IBP

Click or touch the Alarms button on the ECG/IBP tab to show the alarm options for ECG and IBP monitors.

Because the monitor is intended for use with only one animal at a time there are potentially several sources of heart rate data, not least of which are ECG and IBP. When setting a heart rate alarm it is possible to choose to derive the heart rate from either ECG or IBP depending which is in use or which is more reliable in the current circumstances.

The screenshot shows the 'Alarms' settings screen for the ECG/IBP tab. The screen has a top navigation bar with tabs: ECG/IBP, CO2, Oxygen, Agent, N2O, SpO2, AP, Temp, and Traces/Trends. The 'ECG/IBP' tab is selected. On the left, there are two buttons: 'Trace' and 'Alarms', with 'Alarms' being highlighted. Below these buttons is a checkbox labeled 'Synchronise' which is checked. The main area is titled 'Alarm Limits' and contains four sections: 'Heart Rate', 'MAP', 'Systolic', and 'Diastolic'. Each section has 'High' and 'Low' value fields, an 'Enable' checkbox, and a 'Source' selection for Heart Rate (ECG or IBP). The 'Heart Rate' section shows High at 55 bpm and Low at 0 bpm, with ECG selected as the source. The 'MAP' section shows High and Low at 0 mmHg. The 'Systolic' section shows High and Low at 0 mmHg. The 'Diastolic' section shows High and Low at 0 mmHg. A 'Close' button is located at the bottom right.

### Alarm Limits

This group of alarms includes alarms for Heart Rate, Systolic, Diastolic and Mean Arterial Pressure. To set any alarm value, click or touch the associated edit box. The numerical keypad will appear on the right of the screen. Use this or the main keyboard to enter the numerical value and then click or touch Enter.

### Heart Rate

Upper and lower alarm values can be set for the heart rate and these can be independently enabled by touching or clicking the Enable box next to the appropriate Alarm edit box.

### Systolic

Set high and low limits for the systolic blood pressure value. High and Low alarms can be independently enabled.

### Diastolic

Set high and low limits for the diastolic blood pressure value. High and Low alarms can be independently enabled.

### Mean Arterial Pressure

Set high and low limits for the mean arterial blood pressure value. High and Low alarms can be independently enabled.

## USING THE TAFONIUS SOFTWARE

### The CO2 Attribute tab

#### Properties

**Units:** Set the displayed units as Volume %, millimetres of mercury (mmHg) or as kilopascals (kPa). When units are changed, there is an automatic change in the units applied to alarms so that the corresponding value in the new units is used.

**Trace:** Select whether to display the Capnogram as a filled trace or as a simple line curve.

**Scale:** Fixes the vertical scale of the CO2 trace.

**Beep on Breath:** Select this tick box to produce an audible beep with every detected breath.

#### Alarm Options

To set Alarms associated with CO2 measurement, click or touch the Alarms button


Alarm Limits				Enable
Inspired	High	---	%	<input type="checkbox"/>
Expired	High	---	%	<input type="checkbox"/>
	Low	---	%	<input type="checkbox"/>
Breathing Rate	High	---	bpm	<input type="checkbox"/>
	Low	---	bpm	<input type="checkbox"/>

To set an alarm, place a tick in the enable box and then enter the appropriate alarm value.



## USING THE TAFONIUS SOFTWARE

### The Airway Pressure (AP) Attribute tab



#### Sensitivity


This slider sets the scales of the Airway Pressure trace. The value next to the pointer represents the maximum positive height of the trace window.

#### Sweep

This slider sets the sweep speed of the Airway Pressure trace, which is the speed at which the trace is drawn across the screen.

#### Airway Pressure Alarms

To access the Airway Pressure Alarms, click or touch the Alarms button



The CPAP alarm allows an alert condition when the measured CPAP pressure deviates outside of the required CPAP value. Alarm units are in percentage values of the set CPAP value.

## USING THE TAFONIUS SOFTWARE


### The Oxygen Attribute tab

ECG/IBP CO2 **Oxygen** Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Sweep**

Slow  Fast

**Scale**

☐ 0-50%

☒ 0-100%

Close

#### Sweep

Use this slider control to set the sweep speed of the Oxygen trace across the screen.

#### Scale

Use this control to set the full scale height of the Oxygen trace window. Note that there is no shift feature with this trace. To view Oxygen values in excess of 50%, the 0-100% scale range must be used.

#### Alarm Options

Click or touch the Alarms button to show the Alarm options for Oxygen.

ECG/IBP CO2 **Oxygen** Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Alarm Limits**

			Enable
Inspired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
Expired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
FiO2	+/-	<input type="text"/> %	<input type="checkbox"/>

Close

The Oxygen Alarm options allow the alarm limits for Inspired and Expired oxygen to be set. Alarms can be set for High and Low Inspired and Expired values. Click or touch the Enable box to enable the alarm and then enter the alarm value by clicking or touching the appropriate alarm edit box.

**FiO2** - This option is not yet implemented and will be available in the future.

# USING THE TAFONIUS SOFTWARE

## The SPO2 Attribute tab

ECG/IBP CO2 Oxygen Agent N2O **SpO2** AP Temp Traces/Trends

Trace

Alarms

Sweep

Slow Fast

Sound

☐ Beep On SpO2 Pulse

Close

### Sweep

Use the Sweep slider control to change the speed at which the SpO2 trace moves across the screen

### Sound

Click or touch the box next to “*Beep On SpO2 Pulse*” to enable audible beeps associated with each detected pulse waveform. This will over-ride any other beep selection.

### Alarm Options

Click or touch the Alarms button to show the Alarm options for SpO2.

ECG/IBP CO2 Oxygen Agent N2O **SpO2** AP Temp Traces/Trends

Trace

Alarms

Alarm Limits

				Enable
Pulse Rate	High	<input type="text"/>	bpm	<input type="checkbox"/>
	Low	<input type="text"/>	bpm	<input type="checkbox"/>
Saturation	Low	<input type="text"/>	%	<input type="checkbox"/>

Sound

☐ Beep On SpO2 Pulse

Close

### Alarm Limits

Click or touch the enable box to enable the alarm and then enter the appropriate value.

### Sound

Click or touch the tick box to select audible beeps every time a pulse wave on the SPO2 monitor is detected.

## USING THE TAFONIUS SOFTWARE

### The Agent Attribute tab

ECG/IBP CO2 Oxygen **Agent** N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**In Use**

- ☒ No Agent
- ☐ Halothane
- ☐ Enflurane
- ☐ Isoflurane
- ☐ Sevoflurane
- ☐ Desflurane

**Sweep**

Slow  Fast

**Scale**

- ☒ 0-5%
- ☐ 0-10%

Close

**In Use:** The IRMA OR detector used for agent monitoring does not have automatic agent ID so it is necessary to set the agent being used. Select the appropriate radio button next to the agent in use. This change will be reflected in the Agent button on the screen.

### Sweep

Use the Sweep slider control to change the speed at which the Agent trace moves across the screen

### Scale

Use this control to set the full scale height of the Agent trace window. Note that there is no shift feature with this trace. To view Agent values in excess of 5%, the 0-10% scale range must be used.

### Alarm Options

Click or touch the Alarms button to show the Alarm options for Agent.

ECG/IBP CO2 Oxygen **Agent** N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

**Isoflurane Alarm Limits**

			Enable
Inspired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>
Expired	High	<input type="text"/>	<input type="checkbox"/>
	Low	<input type="text"/>	<input type="checkbox"/>

Close

**Alarm Limits:** Set the alarm limits required for a particular agent here. Note that different alarm limits can be set for each and every agent. Click or touch the Enable box to enable the alarm and then enter the appropriate value.

## USING THE TAFONIUS SOFTWARE

### The N2O Attribute tab

ECG/IBP CO2 Oxygen Agent **N2O** SpO2 AP Temp Traces/Trends

Trace

Alarms

**Sweep**

Slow Fast

**Scale**

☐ 0-50%

☒ 0-100%

Close

### Alarm Options

Click or touch the Alarms button to show the Alarm options for N2O.

ECG/IBP CO2 Oxygen Agent **N2O** SpO2 AP Temp Traces/Trends

Trace

Alarms

**Alarm Limits**

Inspired	High	---	%	Enable <input type="checkbox"/>
Expired	High	---	%	Enable <input type="checkbox"/>

Close

### Alarm Limits:

This tab allows the alarm limits to be set. Alarms can be set for High Inspired and Expired values. Click or touch the Enable box to enable the alarm and then enter the alarm value by clicking or touching the appropriate alarm edit box.

## USING THE TAFONIUS SOFTWARE

### The Temperature Attribute tab

ECG/IBP CO2 Oxygen Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

Time

- ☐ 10min
- ☐ 20min
- ☒ 1hr
- ☐ 2hrs
- ☐ 5hrs
- ☐ 10hrs

Temperature Scales

- ☐ Fahrenheit
- ☒ Celcius

Close

#### Time

Selecting an option here sets the amount of data displayed in the trace window. By selecting a long Time period, trends can easily be visualised. By choosing a short time period, details in trace changes are more easily observed.

**Temperature Scales:** Choose to use either Fahrenheit or Celsius as the temperature scale. Alarm values change accordingly.

#### Alarm Options

Click or touch the Alarms button to show the Alarm options for Temperature.

ECG/IBP CO2 Oxygen Agent N2O SpO2 AP Temp Traces/Trends

Trace

Alarms

T1 T2 Delta

Alarm Limits

		Enable
High	<input type="text"/>	<input type="checkbox"/>
Low	<input type="text"/>	<input type="checkbox"/>

Alarm Silence


Close

**Alarm Limits:** There are 3 alarm options: Alarms for T1 temperature probe, Alarms for T2 temperature probe and Alarms for the Delta temperature. (The Delta temperature is the difference in temperature between the two temperature readings of T1 and T2). Click or touch the enable box to enable the alarm and then enter the appropriate value.

## USING THE TAFONIUS SOFTWARE

**Alarm – Silence:** Use this button to silence a temperature alarm. An alarm will sound when a temperature probe is discovered to be absent. When starting the software if no temperature sensors are present an alarm may sound and the border of the temperature gauge will flash red. Touch the temperature gauge to open this Attribute tab and use the Silence button to silence the alarm.

### The Traces/Trends Attribute tab



For a detailed description of how traces are managed and moved around the screen, see the section entitled "Adding and Removing Traces" at the start of the section above.

## **USING THE Tafoneius SOFTWARE**

### **Alarms, Alerts and Critical Error Messages**

DRAFT



# USING THE Tafonius SOFTWARE

## Creating an Anaesthetic Record Chart

The Tafonius software creates an automatic record that can be viewed as an Anaesthetic Record Chart. Every 5 minutes monitoring data is recorded and entered on the chart in a similar manner to a manual record. Information about the patient, surgeon, anaesthetist and other details can all be recorded.


All non-clinical information is entered via the Patient Data button. This process is covered fully in the section entitled "Entering Initial Patient Information"

During a patient Session, there are effectively two recorded sets of data. One set of data is saved as the Slow Data file. This records periodically (default is every 5 seconds) all of the information relating to the running of Tafonius. For a more details explanation of the Slow Data files see the section entitled "Data Recording using the Slow Data File".

The other set of data is data saved as the Anaesthetic Record file. This file consists of data saved every 5 minutes as plotted or tabled data as well as header information such as name of Anaesthetist or Surgeon.

## Altering the Start Time of the Anaesthetic Record chart

It is not uncommon to prepare Tafonius some time in advance of connection to the patient. During this time there will be continual recording of non-meaningful data. To avoid this, the Start time for the Anaesthetic Record data can be set using the Anaesthetic Record Start button. Pressing this button will reset the start time for the Anaesthetic Record data ONLY. Slow data information is not affected. The Anaesthetic Record Start button is found in the Data Logging tab of the Patient Data dialogs. The following dialog will appear:



This procedure can be repeated as many times as required, noting that each time, data prior to selecting this option will no longer be available.

After this creation of the Anaesthetic Record chart is automatic and the chart can be viewed from the Tafonius Shell via the View Data button.

## Viewing the Anaesthetic Record chart

The Anaesthetic Record chart can only be viewed from the Tafonius Shell, whilst Tafonius is **NOT** running. Click or touch the View Data button in the Tafonius Shell. Use the Record Picker to select a patient and then click Select. Then highlight the AR file and click or touch Select. This will open a pdf version of the Anaesthetic Record in Adobe Viewer.

# USING THE TAFONIUS SOFTWARE

## Data Recording using the Slow Data File

The Slow Data file is an automatically generated file that records all information gathered during a Session. Data contained in this file includes all the monitoring information in terms of numerical and trend data. This periodic data (Default is every 5 seconds) also contains all events that occur during use. If an alarm condition occurs this is recorded in the Slow Data file. Similarly a change of anaesthetic agent type or the level of the MWPL setting is recorded in the Slow Data file. This raw file is non-editable so that the user can be sure that all events that occurred are recorded. For a full description of the data saved in the Slow Data file, see Appendix A, Slow Data Information.

The Slow Data file can be exported to a Flash disk using the Export data function accessed from the Tefonius shell. See the section entitled "Export Data" later in this manual.

In addition to automatic events, user events can be added to the Slow Data file. These can be added by clicking or touching the Events button at the bottom right hand side of the Tefonius screen. These events will appear in the Slow Data file as well as being used for the events shown in the Anaesthetic Record chart. For more information on entering manual Events see the section entitled "Adding Events During a Procedure" later in this manual.

## Changing the interval for Periodic Data in the Slow Data File

The interval for periodic data can be changed in the Data Logging tab of the Patient Data dialog. Click or touch the Patient Data button and then select the Data Logging tab. In the Slow Data Logging edit box enter the value in seconds to be used as the logging interval. Note that this value will have an effect on the final file size. Logging every second will create a substantially larger file than every 5 seconds. Note also the effect of the logging interval as explained below.

## Effect of the Logging Interval on the Slow Data entries


Whatever the interval period between logging of entries, automatic event entry will always be logged. For example, if the Ventilate button is pressed an automatic event will be generated and added to a temporary list of events. When the next logging episode occurs, all events stored since the previous logging event will be entered into the SD file. This is different from the physiological monitoring data which is simply grabbed at the instant of the periodic logging.

Care should be used when interpreting periodic logged data because events that appear at the same logging could be separated by up the periodic logging interval.

## USING THE TAFONIUS SOFTWARE

### Adding Events during a procedure

Events can be added to the Slow Data (SD) record during a procedure. Events can be pre-existing in a drop down list or may be added as a new text entry. To add an event, click or touch on the “Events” button at the bottom right hand side of the screen.



This will open a new editor panel at the bottom of the screen.


A screenshot of the 'Event List' editor panel. It has a title bar 'Event List'. Below it is a table with three columns: 'Date/Time', 'Event Description', and 'Status'. The first three rows are populated with: '06/01/2009 15:50 - Fluid Administration:', '06/01/2009 15:49 - Intubation:', and '06/01/2009 15:49 - Anaesthetic Induction:'. To the right of the table are 'Add' and 'Close' buttons.

Click or touch on the “Add” button and then click or touch on the arrow of the Standard Events drop down box. In the example below, Anaesthetic Induction has been chosen from the Drop Down List.

A screenshot of the 'Event Details' editor panel. It has a title bar 'Event Details:'. Below it are two dropdown menus for 'Date' (showing '05 November, 2008') and 'Time' (showing '17:40:58'). Below these is a text input field with 'Anaesthetic Induction: |'. To the right of the text field is an 'Accept' button. At the bottom left is a 'Clear' button. At the bottom center is a 'Standard Events:' label followed by a dropdown menu showing 'Anaesthetic Induction' and an 'Add/Remove Standard Events' button. At the bottom right is a 'Cancel' button.

## USING THE TAFONIUS SOFTWARE


Standard Events can be added or removed from the Drop Down List by using the “Add/Remove Standard Events” button. Clicking or touching this button opens a new editor which lists all the current Standard Events. These can be Edited or Deleted using the Edit/Delete buttons. To add new items to the Standard Event list, click or touch the Add button and then type your new Event title.



The 'Standard Events' editor interface features a title bar 'Standard Events'. On the left, there are two buttons: 'Up' and 'Down'. The central area is a list box containing the following items: 'Start Procedure', 'Fluid Administration', 'Premed', 'Anaesthetic Induction', 'End Anaesthesia', 'Arrival At Recovery Site', 'Extubation', and 'Sternal'. To the right of the list box are four buttons: 'Add', 'Edit', 'Delete', and 'Close'.

Once a new Event has been entered it will appear last on the list. To change the position of this new item in the list, use the UP and Down buttons to move the item. Once the new Event and the order of Events has been set, click or touch the Close button to leave the Standard Events editor.

To add your own entry click or touch in the Event Details box and start typing



The 'Event Details' editor interface has a title bar 'Event Details:'. It includes two date/time pickers: '05 November, 2008' and '17:50:45'. Below these is a large text input field with the placeholder text 'Add a custom event here...'. To the right of the text field is an 'Accept' button. At the bottom left is a 'Clear' button. At the bottom center is a 'Standard Events:' label followed by a dropdown menu. Below the dropdown menu is an 'Add/Remove Standard Events' button. At the bottom right is a 'Cancel' button.

Click or touch “Accept” to add it to the SD file.

Use the “Clear” button to remove an entry from the Event Details edit box.

# USING THE Tafonius SOFTWARE

## The Tafonius Shell

The Tafonius Shell is the starting point for using the Tafonius Ventilator software and also the starting point for various related utilities.



The buttons:


### **Tafonius**

Launches the Tafonius ventilation and monitoring software

### **View Data**

This opens a browser for selecting a patient. All data records are shown in date order with the most recent data at the top of the list. Choose the Patient Record you would like to View and click or touch the Select button. Each patient will have a Slow Data file (SD file) and an Anaesthetic Record file (AR file).

## USING THE TAFONIUS SOFTWARE



To view the Slow Data information, highlight the appropriate file and then touch or click the Select button. The Slow Data is viewed in a spreadsheet format and can be used to show monitoring trends. By touching any of the measured parameter columns (e.g. ECG heart rate, Pulse-Ox saturation) a trend view of the data over the course of the procedure is shown at the bottom of the screen. Touching the column again hides the trend view. For greater versatility, export the data to a USB stick and view the data in a commercial spreadsheet of your choice. All Slow Data files are saved as CSV or comma separated value files.


To view the Anaesthetic Record chart, highlight the desired AR file and then touch or click Select. A new window opens to show the Anaesthetic Record pdf file.

For printing and archiving of either the Slow Data or the Anaesthetic Record use the Export Data button from the TAFONIUS Shell to save the data to a USB flash drive and then open the files in your host computer.

## USING THE TAFONIUS SOFTWARE

### Export Data

This button opens a browser to allow you to choose the data to export. All data in the browser is listed in date order with the most recent files at the top of the browser.




First select the patient that you wish to export data for and then click or touch the Next button. Then choose your destination drive (usually a USB flash drive named D:\) and touch or click the Export button. The Slow Data and the Anaesthetic Record data will be copied to your chosen destination.

## USING THE Tafonius SOFTWARE

### Options

This opens a small window with 3 facilities.



### Touch Screen Alignment

Touching or clicking this button starts a utility to align the touch screen to your computers screen. Simply touch the red circles in the way you would normally use the touch screen and follow the prompts to reset the alignment. If you make a mistake simply wait and the utility will time out and be cancelled. You can then try again.

### Recover Anaesthetic Records

Touching or clicking this button will create a new Anaesthetic Record chart if the original has been destroyed or corrupted for some reason.

### Company Name

Enter your practice or company name here. This will appear on the Tafonius shell as well as at the top of the Anaesthetic Record chart.

### Shutdown

This will shut down the PC and put the PC in an off condition. The PC can be started again by either turning on the Tafonius ventilator or by pressing the large ON/OFF button on the side of the PC.



# USING THE TAFONIUS SOFTWARE

## Appendix A : SLOW DATA INFORMATION

Monitoring data recorded in the Slow Data File includes the following:

### User control settings

CPAP/PEEP	Maximum PEEP or CPAP reading cm H <sub>2</sub> O
TV	Set Tidal Volume L
TVcc	Compliance compensated tidal volume L
RR	Respiratory Rate – breaths per minute
IT	Inspiratory Time in milliseconds
IP	Inspiratory Pause time in milliseconds

### Resultant parameters calculated from user settings

MV	Minute Volume L
I Flow	Inspiratory Flow L/min
I:E Ratio	Inspiratory:Expiratory ratio
Exp Time	Expiratory time in milliseconds

### Measured during Standby

TV	Spontaneous breath Tidal Volume in mls
Insp Time	Inspiratory time in ms
Exp Time	Expiratory time in ms

### Resultant parameters calculated from MEASURED values during Standby

RR	Respiratory rate in breaths per minute
MV	Minute Volume in litres
Iflow	Inspiratory flow rate in L/min
I:E Ratio	Inspiratory:Expiratory ratio

### Airway Pressure Measurements

CPAP/PEEP	CPAP(Standby) or PEEP(Ventilate) values in cm H <sub>2</sub> O
PIP	Peak Inspiratory Pressure in cm H <sub>2</sub> O

### Gas analysis Measurements

Insp CO <sub>2</sub>	Fractional %Volume measurement of inspired CO <sub>2</sub>
Exp CO <sub>2</sub>	Fractional %Volume measurement of expired CO <sub>2</sub>
Resp Rate	Breathing rate in breaths per minute
Insp Agent	Inspired Agent Fractional %Volume measurement
Exp Agent	Expired Agent Fractional %Volume measurement
Insp O <sub>2</sub>	Inspired Oxygen %Volume measurement
Exp O <sub>2</sub>	Expired Oxygen %Volume measurement
Atmospheric Pressure	Current Atmospheric pressure in mm Hg

### Nonin SpO<sub>2</sub> Measurements

SPO <sub>2</sub> HR	Heart rate from pulse-oximeter device
Saturation	Percentage oxygen saturation of arterial blood

## USING THE TAFONIUS SOFTWARE

### IBP Measurements

Mean	Mean arterial blood pressure (mean 4 average)
Systolic	Systolic blood pressure (mean 4 average)
Diastolic	Diastolic blood pressure (mean 4 average)
Heart Rate	Patient heart rate from blood pressure device

### Temperatures

T1	Temperature from Thermometer 1 in either °C or °F
T2	Temperature from Thermometer 2 in either °C or °F

### ECG HR

Heart rate from ECG monitoring device

### Battery Data

Batt1	Voltage level of battery 1 (lower 12v battery)
Current1	Charging rate of battery 1 in amps
Batt2	Voltage level of battery 2 (upper 12v battery)
Current2	Charging rate of battery 2 in amps

### Piston Position

Units	Shows the piston position in absolute units
-------	---

In addition to these headings in the file, there will be additional events recorded. The table below shows some typical events that have been recorded.

[08:09:33.21 - Preset Loaded: Custom Preset]  
[08:09:33.191 - Leak 600.0, Compliance Low 27.4, Compliance High 63.3]  
[08:09:33.201 - MWPL: 35]  
[08:09:33.311 - Cylinder Open]  
[08:09:33.702 - FMSO ON: 65535ms]  
[08:09:33.722 - Init Requested]  
[08:09:35.695 - Auxiliary Controller - No Watchdog]  
[08:09:35.705 - Auxiliary Controller - No Position data]  
[08:09:35.895 - Auxiliary Controller - Watchdog data established]  
[08:09:37.237 - FMSO ON: 65535ms]  
[08:09:37.247 - Init Requested]  
[08:09:37.778 - Auxiliary Controller - Position data established]  
[Thursday, March 27, 2008 08:09:38] Records Started - Tafarius Version 0.1.1092.0.

Every entry has a time associated with it and describes the event that occurred.